

NDiff[™] NEURO-2 MEDIUM SUPPLEMENT (200X)

CATALOG NUMBER:	SCM012	QUANTITY:	5 mL
LOT NUMBER:		CONCENTRATION:	200X
DESCRIPTION:	NDiff [™] Neuro-2 Medium Supplement is a serum-free supplement developed for the <i>in vitro</i> differentiation of murine embryonic stem (ES) cells into post-mitotic neurons particularly via monolayer differentiation ¹⁻³ . This product may also be used in the derivation, propagation and maintenance of mouse NS neural stem cells ⁴⁻⁶ .		
FORMULATION:	NDiff [™] Neuro-2 Medium Supplement contains Insulin, Apo-Transferrin, Sodium Selenite, Putrescine, and Progesterone. This media supplement is supplied as a 200X stock and should be added to the cell culture medium at a final concentration of 1X.		
FORMAT:	Frozen liquid.		
STORAGE/HANDLING:	Maintain at -20 °C until expiration date. Protect from light. Once thawed, aliquot unused portion into smaller volumes and store at -20 °C until future use. Once added to the cell culture medium, the product is stable at 4 °C for four weeks. Product is supplied non-sterile. Complete medium must be filtered prior to use.		
REFERENCES:	<ol style="list-style-type: none">1. Nichols J, Ying QL. (2006) Derivation and propagation of embryonic stem cells in serum- and feeder-free culture. <i>Methods Mol Biol.</i> 329:91-8.2. Ying QL, Smith AG. (2003) Defined conditions for neural commitment and differentiation. <i>Methods Enzymol.</i> 365:327-41.3. Ying QL, Stavridis M, Griffiths D, Li M, Smith A. (2003) Conversion of embryonic stem cells into neuroectodermal precursors in adherent monoculture. <i>Nat Biotechnol.</i> 21(2):183-6.4. Conti L, Reitano E, Cattaneo E. (2006) Neural stem cell systems: diversities and properties after transplantation in animal models of diseases. <i>Brain Pathol.</i> 16(2):143-54.5. Pollard SM, Conti L, Sun Y, Goffredo D, Smith A. (2006) Adherent neural stem (NS) cells from fetal and adult forebrain. <i>Cereb Cortex.</i> 16 Suppl 1:i112-20.6. Conti L, Pollard SM, Gorba T, Reitano E, Toselli M, Biella G, Sun Y, Sanzone S, Ying QL, Cattaneo E, Smith A. (2005) Niche-independent symmetrical self-renewal of a mammalian tissue stem cell. <i>PLoS Biol.</i> 3(9):e283.		

© 2006: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.