



- ▶ High recovery
- ▶ High throughput
- ▶ Wide process flexibility

# ProSep<sup>®</sup>-PB Chromatography Media

***For the affinity chromatography of a variety of glycoproteins, nucleic acids, and other biomolecules***

ProSep-PB media is designed for the affinity chromatography of a variety of glycoproteins, nucleic acids, and other biomolecules containing 1,2 diol functionality. ProSep-PB media can be used in the purification of glycoproteins, nucleotides and viruses.

ProSep-PB media is composed of a synthetic m-aminophenyl ligand immobilized on controlled pore size glass beads. The non-compressible nature of the media and the low back pressure generated, allows efficient media utilization at all scales.

Phenyl boronate ligand mediated chromatography differs in one fundamental way from most other ligands used for affinity chromatography. Whereas most ligands bind to a particular binding site on a protein by a mixture of non-covalent interactions, phenyl boronate interacts predominantly by forming a temporary covalent bond with 1,2-cis-diol groups. The boronate ligand will bind to any molecule containing the appropriate group, for example, erythropoietin (EPO), which is highly glycosylated. Steric configuration is only important in how it affects the accessibility of the 1,2-cis-diol groups. The generality of the binding mechanism makes ProSep-PB a versatile chromatography medium.

# Examples of Affinity Chromatography Using Phenyl Boronate as a Ligand

## TARGET MOLECULE

<b>Carbohydrates</b>	Arabinose <sup>7</sup> , fructose <sup>1,7,20</sup> , galactose <sup>4,20</sup> , glucose <sup>1,4,7,20</sup> , inositol <sup>1,20</sup> , lactose <sup>4,7</sup> , mannitol <sup>1,20</sup> , mannose <sup>1,7,14</sup> , ribose <sup>1,4,7</sup> , sorbitol <sup>1,20</sup> , sucrose <sup>1,20</sup>
<b>Nucleosides/Nucleotides</b>	Adenosine <sup>1,2,3,7,11,12,20</sup> , deoxyadenosine <sup>1,7,16,20</sup> , 5'-AMP <sup>2,6,7,20,22</sup> , 5'-ADP <sup>2,6,7,20,22</sup> , 5'-ATP <sup>2,6,7,20,22</sup> , cytidine <sup>1,2,3,11,20</sup> , guanosine <sup>1,3,11,20</sup> , thymidine <sup>1,2,20</sup> , uridine <sup>1,2,3,11,20</sup> , NAD(H) <sup>15,22</sup> , NADP(H) <sup>20,22</sup>
<b>Nucleic Acids</b>	Oligonucleotides <sup>2,5,18</sup> , tRNA aminoacylated from nonaminoacylated forms <sup>5,8,18,19</sup>
<b>Peptides and Proteins</b>	Erythropoiten (EPO) <sup>23</sup> , membrane glycoproteins <sup>21</sup>
<b>Other Small Molecules</b>	Dopamine <sup>10,13,17</sup> , 3-methoxytyramine <sup>13</sup> , norepinephrine <sup>10,17</sup> , tyrosine <sup>9,13</sup>



## Specifications

<b>Matrix</b>	Porous glass with particle size in the range $\geq 74 \leq 125$ microns
<b>Working pH Range</b>	pH 2.5 – 8.5
<b>Leachables</b>	0.3ppm Si at pH7, decreasing at lower pH (0.01 ppm at pH 1.5). ProSep-PB media is stable between pH1 and pH9.
<b>Pressure</b>	Incompressible matrix with linear pressure/flow rate characteristics. Maximum operating pressure > 3000 psi (200 bar)
<b>Shelf Life</b>	1 year
<b>Life</b>	ProSep-PB media is stable over repeated operational cycles provided that proper cleaning protocols are observed.

## References

1. Weith, H., Wiebers, J., and Gilham, P. 1970. Synthesis of cellulose derivatives containing the dihydroxyboryl group and a study of their capacity to form specific complexes with sugars and nucleic acid components. *Biochem.* 9:4396
2. Rosenberg, M., Weibers, J., and Gilham, P. 1972. Studies on the interactions of nucleotides, polynucleotides and nucleic acids with dihydroxyboryl-substituted celluloses. *Biochem.* 11:3623
3. Schott, H., 1972. New dihydroxyboryl-substituted polymers for column chromatographic separation of ribonucleoside-deoxyribonucleoside mixtures. *Agnew. Chem.* 84:819
4. Reske, K., and Schott, H. 1973. Column chromatographic separation of neutral sugars on a dihydroxyboryl-substituted polymer. *Agnew. Chem.*, 85:417
5. Schott, H., et al. 1973. A dihydroxyboryl-substituted methacrylic polymer for the column chromatographic separation of mononucleotides, oligonucleotides, and transfer ribonucleic acid. *Biochem.*, 12:932
6. Moore, E., et al. 1974. Separation of ribonucleotides and deoxyribonucleotides on columns of borate covalently linked to cellulose. Application to the assay of ribonucleoside diphosphate reductase. *Biochem.*, 13:2904
7. Yurkevich, A., et al. 1975. Study of the interaction of polyols with polymers containing N-substituted [(4-boronophenyl) methyl] - ammonio groups. *Carbohyd. Res.*, 43:215
8. McCutchan, T., Gilham, P., and Soll, D. 1975. An improved method for the purification of RNA by chromatography on dihydroxyboryl substituted cellulose. *Nucleic Acids Res.*, 2:853.
9. Ellinger, C., Chan, B., and Stanley, W. 1975. P-Vinylbenzeneboronic acid polymers for the separation of vicinal diols. *J. Chrom.*, 104:57
10. Higa, S., et al. 1977. Isolation of catecholamines in biological fluids by boric acid gel. *Anal. Biochem.*, 77:18
11. Uziel, M., Smith, L., and Taylor, S. 1976. Modified nucleosides in urine: selective removal and analysis. *Clin. Chem.*, 22:1451
12. Davis, G., et al. 1977. HPLC separation and quantitation of nucleosides in urine and some other biological fluids. *Clin. Chem.*, 23:1427
13. Hansson, C., et al. 1978. Chromatographic separation of catecholic amino acids and catecholamines on immobilised phenylboronic acid. *J. Chrom.*, 161:352
14. Wulff, G., and Vesper, W. 1978. Preparation of chromatographic sorbents with chiral cavities for racemic resolution. *J. Chrom.*, 167:171
15. Maestas, R., et al. 1980. Polyacrylamide-boronate beads saturated with biomolecules: A new general support for affinity chromatography of enzymes. *J. Chrom.*, 189:225
16. Zimmerman, T., Wolberg, G., Duncan, G., and Elion, G. 1980. Adenosine analogues as substrates and inhibitors of S-adenosylhomocysteine hydrolase in intact lymphocytes. *Biochem.*, 19:2252.
17. Gelijkens, C., and Deleenheer, A. 1980. Simple method for isolating free urinary catecholamines by boric acid gel chromatography. *J. Chrom.*, 183:78
18. Pace, B., and Pace, N. 1980. The chromatography of RNA and oligoribonucleotides on borate-substituted agarose and polyacrylamide. *Anal. Biochem.*, 107:128
19. Singhal, R., Bajaj, R., Buess, C., Smoll, D. and Vakharia, V. 1980. Reversed-phase boronate chromatography for the separation of o-methyl-ribose nucleosides and amino-t-RNA's. *Anal. Biochem.*, 109:1
20. Glad, M., et al. 1980. High performance liquid affinity chromatography of nucleosides, nucleotides and carbohydrates with boronic acid substituted microparticulate silica. *J. Chrom.*, 200:254
21. Williams, G., Johnstone, A., Bouriotis, V., and Dean, P. 1981. Affinity chromatography of membrane proteins on dihydroxyboryl-matrix gel. *Biochem. Soc. Trans.*, 9:137
22. Bouriotis, V., Galpin, I., and Dean, P. 1981. Applications of immobilised phenylboronic acids as supports for group specific ligands in the affinity chromatography of dehydrogenases. *J. Chrom.*
23. D. Zanette, A. Soffientini, C. Sottani, E. Sarubbi. 2003. Evaluation of phenylboronate agarose for industrial scale purification of erythropoietin (EPO) from mammalian cultures. *J. Biotech* 101, 275-287

## Ordering Information

Description	Pack Size	Former CAE Reference No.	Catalogue No.
ProSep-PB Media	10 mL	CAE 004124	113 247 324
	50 mL	CAE 004126	113 247 326
	100 mL	CAE 004127	113 247 327
	500 mL	CAE 004129	113 247 329
	1000 mL	CAE 004130	113 247 330
	5 L	CAE 004133	113 247 333
	10 L	CAE 004135	113 247 335

### Discover the More in Millipore™

In every application, every step and every scale, count on Millipore to be everywhere for you — from monoclonals to vaccines, from clinical through pilot to full-scale manufacturing. Our technologies are used by most of the world's major biopharmaceutical companies. But we deliver more than advanced separation, purification, sterilization and quality control products. With Millipore, you get services to optimize and validate your processes, comprehensive resources to streamline and enhance your operation, unmatched know how forged from nearly 50 years' experience — and solutions that integrate it all. For higher yields, improved process economics and faster speed to market, discover the more in Millipore.

### To Place an Order or Receive Technical Assistance

For additional information call your nearest Millipore office:

In the U.S. and Canada, call toll-free  
**1-800-MILLIPORE (1-800-645-5476)**

In the U.S., Canada and Puerto Rico,  
fax orders to **1-800-MILLIFX  
(1-800-645-5439)**

Outside of North America, contact  
your local office. To find the office  
nearest you:

[www.millipore.com/offices](http://www.millipore.com/offices)

Internet: [www.millipore.com](http://www.millipore.com)

Tech Service:

[www.millipore.com/techservice](http://www.millipore.com/techservice)

Millipore and PROSEP are registered trademarks of  
Millipore Corporation.

Lit. No. DS1055EN00 Rev. - 9/03 03-104

© 2003 Millipore Corporation, Billerica, MA  
01821 U.S.A. All rights reserved.

# MILLIPORE