



**NF $\kappa$ B p50/p65**  
(Human/Mouse/Rat)  
**Transcription Factor Assay**  
**Colorimetric**

**Cat. No. SGT510 NF $\kappa$ B p50/p65**

**Cat. No. SGT515 NF $\kappa$ B p50**

**Cat. No. SGT520 NF $\kappa$ B p65**

**FOR RESEARCH USE ONLY**  
**Not for use in diagnostic procedures**

USA & Canada

Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233

Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered Worldwide

[www.chemicon.com](http://www.chemicon.com) • [custserv@chemicon.com](mailto:custserv@chemicon.com) • [techserv@chemicon.com](mailto:techserv@chemicon.com)

This page left blank intentionally

---

## Introduction

The transcription factor NF $\kappa$ B (Nuclear Factor kappa B) is involved in the expression and regulation of a number of important cellular and physiological processes such as growth, development, apoptosis, immune and inflammatory response, and activation of various viral promoters including human immunodeficiency virus long terminal repeats<sup>1,2</sup>. NF $\kappa$ B represents a group of structurally related and evolutionarily conserved proteins related to the proto-oncogene c-Rel with five members in mammals that include Rel (cRel), RelA (p65), RelB, NF $\kappa$ B1 (p50 and its precursor p105), and NF $\kappa$ B2 (p52 and its precursor p100)<sup>1, 2</sup>. NF $\kappa$ B/Rel proteins exist as homo- or heterodimers to form transcriptionally competent or repressive complexes. Although most NF $\kappa$ B dimers are activators of transcription, the p50/50 and p52/52 homodimers can repress the transcription of their target genes<sup>2, 3</sup>. The p50/p65 heterodimer of NF $\kappa$ B is the most abundant in cells<sup>3</sup>.

A critical component in NF $\kappa$ B regulation is the I $\kappa$ B Kinase (IKK) complex<sup>3</sup>. In a majority of unstimulated cells, the NF $\kappa$ B transcription factors exist in their inactive form and are retained in the cytoplasm by the bound inhibitory I $\kappa$ B proteins<sup>4,5</sup>. Upon stimulation by multiple inducers including viruses or cytokines, such as TNF $\alpha$ , IL-1, or PMA, I $\kappa$ B $\alpha$  is rapidly phosphorylated and degraded, resulting in the release of the NF $\kappa$ B complex, most commonly the p105/p65 heterodimer. The p105 subunit is cleaved into its active p50 form. This cleavage exposes the NLS sequence on the p50 subunit. The p50/p65 heteroduplex then translocates to the nucleus where it activates gene transcription. NF $\kappa$ B induces the transcription of its own inhibitor, I $\kappa$ B $\alpha$ , causing an autoregulatory mechanism of NF $\kappa$ B activity and generating the inactive form of NF $\kappa$ B. The newly formed nuclear NF $\kappa$ B-I $\kappa$ B $\alpha$  complexes are then exported out to the cytoplasm, thereby reestablishing the cytoplasmic pool of inactive NF $\kappa$ B complexes primed for another round of activation to take place<sup>2</sup>. The wide variety of genes regulated by NF $\kappa$ B includes those encoding cytokines, chemokines, adhesion molecules, acute phase proteins, and inducible effector enzymes<sup>3</sup>.

---

## Application

The CHEMICON<sup>®</sup> Non-Radioactive NFκB p50/p65 Transcription Factor Assays provides a fast, sensitive method to detect specific transcription factor DNA binding activity in nuclear extracts. This assay combines the principle of the electrophoretic mobility shift assay (EMSA) with the 96-well based enzyme-linked immunosorbent assay (ELISA), enabling manual or high-throughput sample analysis with significantly greater sensitivity than conventional EMSA assays. The entire assay takes less than 4 hours to complete with minimal hands-on time. The versatile set up allows for a flexible assay design and because the binding reaction occurs in solution, various parameters can be optimized such as probe titration, competitive oligonucleotide concentration, or treatment conditions.

---

## Test Principle

The Non-Radioactive NFκB p50/p65 Transcription Factor Assay kit is provided in a 96-well format. During the assay, the Capture Probe, a double stranded biotinylated oligonucleotide containing the flanked DNA binding consensus sequence for NFκB (5'-GGGACTTTCC-3'), is mixed with cellular (nuclear) extract in the Transcription Factor Assay Buffer provided. When incubated together, the active form of NFκB contained in the nuclear extract binds to its consensus sequence. The extract/probe/buffer mixture is then directly transferred to the streptavidin-coated plate. The active NFκB protein is immobilized on the biotinylated double stranded oligonucleotide capture probe bound to the streptavidin plate well, and any inactive, unbound material is washed away. The bound NFκB transcription factor subunits, p50 and/or p65, are detected with specific primary antibodies, a Rabbit anti-NFκB p50 and a Rabbit anti-NFκB p65. A highly sensitive HRP-conjugated secondary antibody is then used for detection. This provides sensitive colorimetric detection that can be read in a spectrophotometric plate reader. Included in the kits are positive cell extract, a non-specific double stranded oligonucleotide, and a specific competitor double stranded oligonucleotide.

*For Research Use Only; Not for use in diagnostic procedures*

*Use of this assay in NFκB-related drug discovery may be covered under U.S. Patent No. 6,150,090 and require a license from Ariad Pharmaceuticals (Cambridge, MA, USA).*

---

## Kit Components

1. Transcription Factor Assay Plate (Part No. 90268): One Streptavidin coated 96-well strip plate (12 removable strips of 8 wells each), sealed in a foil pouch with desiccant.
2. TFA Buffer, 5x (Part No. 90272): One bottle containing 80 mL of 5x Transcription Factor Assay Buffer containing sonicated salmon sperm DNA to block non-specific DNA binding activity.
3. Blocking Reagent (Part No. 90273): One vial containing 12 grams of Blocking Reagent.
4. Primary Antibody (Part No. 90274): (SGT510 and SGT520) One vial containing 15  $\mu$ L Rabbit anti-NF $\kappa$ B p65. Use at 1:1000 dilution in Enhanced Transcription Factor Assay Buffer (see Assay Protocol).
5. Primary Antibody (Part No. 90275): (SGT510 and SGT515) One vial containing 15  $\mu$ L of anti-NF $\kappa$ B p50. Use at 1:1000 dilution in Enhanced Transcription Factor Assay Buffer (see Assay Protocol).
6. Gt X Rb, HRP Secondary Antibody (Part No. 90276): One vial containing 30 $\mu$ L Goat anti-Rabbit, HRP conjugated secondary detection antibody. Use at 1:500 dilution in Enhanced Transcription Factor Assay Buffer.
7. NF $\kappa$ B Capture Probe (Part No. 90269): One vial containing 220  $\mu$ L biotinylated double stranded oligonucleotide probe (1 pmol/ $\mu$ L) containing the wild-type consensus sequence (5'-GGGACTTCC-3') for NF $\kappa$ B.
8. NF $\kappa$ B Specific Competitor Oligonucleotide (Part No. 90270): One vial containing 20  $\mu$ L (10 pmol/ $\mu$ L) specific competitor oligonucleotide. The NF $\kappa$ B Competitor Oligonucleotide has the same consensus sequence (5'-GGGACTTCC-3') as the NF $\kappa$ B Capture Probe but is not biotinylated and will not bind to the plate. It is provided as a specific competitor for transcription factor binding and can be used to demonstrate specificity. This control does not need to be run with each assay performed.
9. TFA Negative Control Probe (Part No. 90271): One vial containing 20  $\mu$ L biotinylated double stranded oligonucleotide probe (2 pmol/ $\mu$ L). Provided as an internal negative control, this non-specific competitor oligonucleotide does not contain the NF $\kappa$ B consensus sequence.
10. HeLa Whole Cell Extract, TNF $\alpha$  Treated (Part No. 90350): One vial containing 25  $\mu$ L TNF $\alpha$  stimulated HeLa whole cell extract. This extract is provided as a positive control and is not intended for extract-to-extract, nor plate-to-plate comparisons. When using this control, you will notice a decrease in signal with repeated freeze/thaw cycles. This component must be stored at -80°C.

11. Substrate (TMB/E) (Part No. 90348): One 12mL bottle (Ready to use).
12. Stop Solution (Part No. 60193): One 12mL bottle containing 0.5M HCl.  
*(Caution: Caustic material. Eye, hand, face, and clothing protection should be worn when handling this material.)*

---

### **Materials Not Supplied**

1. Nuclear Extraction Kit (Catalog No. 2900) or suitable extraction reagents
2. Pipettor (both single and multi-channel)
3. Pipette tips
4. Distilled water
5. Microplate reader (wavelength 450nm/650nm) (650nm optional)
6. Rotator or Plate shaker (optional)
7. Plate cover (optional)

---

### **Storage**

The Transcription Factor Assay is a triple temperature storage assay kit.

- Store kit components at the temperatures indicated on the labels and this insert until their expiration date.
- Store the anti-NF $\kappa$ B p50 and p65 primary antibodies, the secondary anti-rabbit-HRP antibody, the NF $\kappa$ B Capture Probe, the NF $\kappa$ B Competitor Oligonucleotide, and the Negative Control Probe at -20°C.
- Store the Transcription Factor Assay (TFA) Buffer at -20°C. Avoid freeze/thawing of the TFA Buffer. Buffer can be thawed in warm water prior to use. Buffer can be aliquotted and frozen.
- Store the cell lysate at -80°C. Refrain from repeated freeze/thawing. Always keep on ice when using. Stable at -80°C for 6 months.
- Store all other components at 4°C.

---

## Sample Preparation

(Performed Prior to Assay)

### Cytosolic and Nuclear Extraction Procedure

Prior to performing the Chemicon Transcription Factor Assay, nuclear or whole cell protein extraction needs to be performed. To perform a nuclear extraction, we recommend using Chemicon's Nuclear Extraction Kit (Cat. No. 2900), or for a quick procedure, Chemicon offers a Whole Cell Extraction kit (Cat. No. 2910).

If the above referenced kits are not used to perform the protein extraction, the following procedure is recommended. Various other methods could be utilized to prepare cytosolic and nuclear extracts, though the following recommended protocol has been tested and works well with the Chemicon Transcription Factor Assay.

#### Materials Required for Nuclear Extraction:

1. PBS (pH 7.5)
2. Protease Inhibitor Cocktail (in DMSO)
3. Syringe with No. 27 needles (small gauge)
4. Materials for buffers: HEPES, MgCl<sub>2</sub>, KCl, EDTA, NaCl, Triton X-100, DTT, glycerol, and Igepal CA-630 (NP-40 can be used in place of Igepal CA-630).

Make Buffer A and Buffer B required for the nuclear extraction in advance. Store both buffers on ice.

*Note: Fresh DTT and protease inhibitors should be added just prior to use.*

#### Buffer A (Hypotonic Lysis Buffer)

10 mM HEPES (pH 7.9)  
1.5 mM MgCl<sub>2</sub>  
10 mM KCl  
0.5-5 mM DTT  
0.1% Triton X-100  
Protease inhibitor cocktail

*Note: Buffer A can be made at a 10x concentration and stored at 2-8°C without the DTT and Protease inhibitors added.*

*Note: Increased amounts of DTT may be necessary dependent upon cell type and/or treatment.*

**Buffer B (Extraction Buffer)**

20 mM HEPES (pH 7.9)  
1.5 mM MgCl<sub>2</sub>  
0.42 M NaCl  
0.2 mM EDTA  
0.5-5 mM DTT  
1.0% Igepal CA-630  
25% (v/v) glycerol  
Protease inhibitor cocktail

*Note: All work done after cell trypsinization/detachment needs to be performed on ice and/or with chilled buffers. It is imperative that the cell pellets and suspension remain as cold as possible without freezing during the extraction process.*

*Note: Increased amounts of DTT may be necessary dependent upon cell type and/or treatment.*

**A. Cell Culture**

1. Grow cells to desired confluency and add appropriate activators and/or inhibitors.

**B. Cell Disruption**

1. For adherent cells, wash the cells with HBSS (Hank's Balanced Salt Solution) or PBS and then add warmed trypsin to the culture flask(s). Let the trypsin sit for about 2 minutes and shake the cells off. Alternately, cell lifters may be used instead of trypsin. Collect cells and transfer them to a clean centrifuge tube, rinsing the culture flask with ice cold PBS. Centrifuge the sample at 250 x g for 5 minutes. Discard the supernatant and resuspend the cell pellet in 40 mL of ice-cold PBS. Centrifuge the suspension at 250 x g as before. Repeat. Pour off supernatant.
2. Estimate the approximate volume of the centrifuged cell pellet. This value will be needed for determining the amount of the various buffers to add for nuclear extraction.
3. Add 5 cell pellet volumes of ice cold Buffer A (lysis buffer).
4. Resuspend the cell pellet by gently inverting the tube avoiding foam production. Do not vortex.
5. Incubate the cell suspension in ice cold Buffer A for 15 min on ice.
6. Centrifuge the cell suspension for 5 minutes at 250 x g. Discard supernatant and resuspend the cell pellet in two volumes of ice cold Buffer A.

### **C. Cell Lysis**

*Note: All subsequent steps must be performed on ice.*

1. Using a syringe with a small gauge needle (a No. 27 gauge needle works well) slowly draw the cell suspension prepared in Section B, Step 6 into the syringe and then eject the contents with a single stroke. Repeat approximately 5 times (drawing and ejecting). If the cells “clump” and you are not able to draw them into the syringe with a No.27 gauge needle, more non-ionic detergent (Igepal CA-630 or NP40) may need to be added.
2. Centrifuge the disrupted cell suspension at 8,000 x g for 20 minutes.
3. The supernatant contains the cytosolic portion of the cell lysate. Transfer the supernatant to a fresh tube. To keep the cytosolic fraction, snap-freeze and store in aliquots at  $-70^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.
4. The remaining portion (pellet) contains the nuclear portion of the cell lysate.

### **D. Nuclear Extraction**

1. Resuspend the nuclear pellet in 2/3 of the original cell pellet volume (determined in Step B.2) ice cold Buffer B.
2. Use a fresh syringe (with a No. 27 gauge needle), repeat Step C.1. to disrupt the nuclei.

*Note: The nuclear extract sample can be stored at  $-80^{\circ}\text{C}$  at this point if needed.*

3. Use a rotator (low speed) to gently agitate the nuclear suspension at  $4^{\circ}\text{C}$  for 30-60 minutes.
4. Centrifuge the nuclear suspension at 16,000 x g for 5 min.
5. Transfer the supernatant to a fresh tube. This fraction is the nuclear extract.
6. Determine protein concentration and adjust the final protein concentration to 2.5 to 5 mg/mL with ice cold PBS.
7. Snap-freeze the nuclear extract in aliquots and store at  $-80^{\circ}\text{C}$ . Avoid repeated freezing and thawing of nuclear extract.

---

## Assay Protocol

### A. Prepare 1x Transcription Factor Assay Buffer (TFA Buffer)

1. Calculate the total volume of 1x TFA Buffer required: Each well requires 2.7 mL of 1x TFA Buffer (i.e. 8 wells require 21.6 mL of 1x TFA Buffer).
2. Dilute 1 part 5x TFA Buffer into 4 parts distilled water to obtain the desired amount of buffer.

### B. Prepare Enhanced Transcription Factor Assay Buffer (Enhanced TFA Buffer)

1. Take a portion of the 1x TFA Buffer prepared in Step A and use to prepare 1x Enhanced TFA Buffer by adding blocking reagent.
2. Calculate the total volume of Enhanced TFA Buffer required for the assay run: Each well requires 1.5 mL of 1x TFA Buffer supplemented with 0.09 g Blocking Reagent (i.e. 8 wells require 12 mL of enhanced TFA buffer made with 0.72 g of blocking reagent). The Enhanced TFA Buffer will be used for all steps of the assay with the exception of the final wash following secondary antibody incubation. For this wash, use the remaining 1x TFA Buffer prepared in Step A (approximately 1.2 mL per assay/well).
3. Mix the Blocking Reagent with 1x TFA Buffer making sure that the blocking reagent is completely dissolved in the TFA buffer prior to starting the assay.

### C. Sample Incubation

1. Prepare the test sample, negative control sample, positive control sample, and/or competitive oligonucleotide control sample in a 0.5 mL microcentrifuge tube (1.5 mL microcentrifuge tube or 96-well plate for multiple tests) as described in the table on the next page.

Alternatively, the samples may be prepared directly in the TFA plate when added in sequential order as follows:

- 1st.* 1x Enhanced TFA Buffer
- 2nd.* Oligonucleotide probes
- 3rd.* Cell or nuclear extract

**Note:** *Sample extract must always be added to the binding reaction mixture last*

**Note:** *Never use excessive heat to thaw the capture probes or competitor oligonucleotide. Use of heat can destroy product.*

	Enhanced Transcription Factor Assay Buffer (1x)	NFκB Capture Probe (90269)	NFκB Competitor Oligonucleotide (90270)	TFA Negative Control Probe (90271)	Nuclear Extract	HeLa Whole Cell Extract (TNFα Treated) (90350)	Total Volume
<b>Transcription Factor Assay (Normal)</b>	43μL	2μL	–	–	5μL	–	50μL
<b>Transcription Factor Assay— Positive Control</b>	43μL	2μL	–	–	–	5μL	50μL
<b>Transcription Factor Assay— Specific Competitor Control</b>	41μL	2μL	2μL	–	5μL	–	50μL
<b>Transcription Factor Assay— Negative Control</b>	44μL	–	–	1μL	5μL	–	50μL

**a. Positive Control:**

The positive control assay utilizes TNFα treated HeLa whole cell extract and should be performed to ensure that the assay is performing correctly. Under normal conditions, this assay will result in signal from the sample well(s). It is not necessary to add the two control probes (the NFκB Competitor Oligonucleotide and the TFA Negative Control Probe) to the positive control extract, unless desired.

**b. Competitor Control:**

The specific NFκB Competitor Control oligonucleotide is provided to ensure that the NFκB complex is binding the probe DNA in a sequence specific manner. It is an unlabeled competitor oligonucleotide containing the identical NFκB consensus sequence as the capture probe. This control will compete with the capture probe for NFκB binding. With a typical nuclear extract, this assay setup should greatly diminish the signal intensity. However, assays may be performed using varying amounts of the NFκB capture probe and/or specific competitor oligonucleotide as desired. Alternatively, the investigator may choose to incorporate other competitor oligonucleotides into the assay. *Keep in mind that any increase in volume of added competitor oligonucleotide must be compensated for by an adjustment in the amount of 1x TFA Buffer in the assay setup such that the final volume equals 50μL.*

**c. Negative Control:**

The TFA Negative Control Probe is used without the addition of the NFκB Capture Probe. The TFA Negative Control Probe is used to ensure that the signal obtained from the normal assay setup is specific to NFκB binding to its consensus sequence. Any signal obtained with this assay setup will reflect non-specific binding of NFκB to the capture probe.

2. For samples not prepared directly in the TFA plate well, transfer the entire 50 μL reaction to an individual well on the TFA plate. *Omit this step for samples prepared directly on the TFA plate.* Place any remaining unused assay strips from the plate at 4°C for future assay use.
3. Incubate in plate well for 1-2 hours at room temperature. A 1 hour incubation is sufficient in most cases, but a 2 hour incubation can result in a slightly higher overall signal as well as a slightly higher signal to noise ratio, especially when using less sample. Cover the sample wells during this and all subsequent non-wash incubations.
4. Wash wells three times with 150 μL of the 1x Enhanced TFA Buffer. Incubate each wash as follows: 1<sup>st</sup> wash, 30 seconds; 2<sup>nd</sup> wash, 1 minute; 3<sup>rd</sup> wash, 3 minutes. Following each wash, invert the plate over a basin to empty the wells and gently tap the inverted plate a few times onto an absorbent pad. Completely remove the buffer from the wells following the final wash.

*Note: Mild agitation or rotating during this and all subsequent washes may give cleaner results, but is not necessary.*

**D. Primary Antibody Incubation**

1. Dilute primary antibody in 1x Enhanced Transcription Factor Assay Buffer.
  - a. Prepare enough for 100 μL of diluted antibody per well tested.
  - b. Separately dilute the anti-NFκB p50 primary antibody 1:1000 and/or the anti-NFκB p65 primary antibody 1:1000 in 1x Enhanced TFA buffer.
2. Add 100 μL of diluted primary antibody to each assay well. Incubate for 60 minutes at room temperature.

3. Wash wells three times with 150  $\mu\text{L}$  of the 1x Enhanced Transcription Factor Assay Buffer (as in Step C.4.). Incubate each wash as follows: 1<sup>st</sup> wash, 30 seconds; 2<sup>nd</sup> wash, 1 minute; 3<sup>rd</sup> wash, 3 minutes. Completely remove the buffer from the wells following the final wash.

#### **E. Secondary Antibody Incubation**

1. Dilute the anti-rabbit IgG-HRP conjugated secondary antibody in 1x Enhanced Transcription Factor Assay Buffer at a 1:500 dilution.
  - a. Prepare enough for 100  $\mu\text{L}$  of diluted secondary antibody per assay well tested.
2. Add 100  $\mu\text{L}$  of the diluted secondary antibody to each well and incubate at room temperature for 30 minutes.

*Note: At this time the TMB/E should be removed from 2° to 8°C and allowed to equilibrate to room temperature during this 30 minute incubation. For increased sensitivity warm to 37°C prior to use.*

3. Wash wells four times with 250  $\mu\text{L}$  of 1x Transcription Factor Assay Buffer (**WITHOUT BLOCKING REAGENT**) per well. Incubate each wash as follows: 1<sup>st</sup> wash, 30 seconds; 2<sup>nd</sup> wash, 2 minutes; 3<sup>rd</sup> wash, 3 minutes; 4<sup>th</sup> wash, 4 minutes.

*Note: Using 1x Enhanced TFA Buffer, which contains Blocking Reagent, in this step will increase background signal.*

#### **F. Color Development**

1. Add 100  $\mu\text{L}$  of the pre-equilibrated TMB/E substrate to each assay well being tested. Cover the plate and incubate at room temperature (or 37°C for increased sensitivity) for approximately 10 minutes. The development time may vary according to laboratory conditions. Therefore, we recommend that you closely monitor the speed at which color is generated during the 10 minute incubation. The assay may need to be stopped prior to 10 minutes to prevent signal saturation.


*Note: The development time may be extended to as much as 30 minutes for very low detection levels.*

2. Add 100  $\mu\text{L}$  of the Stop Solution to each well.
3. Measure the absorbance of the samples using a spectrophotometric microplate reader set at  $\lambda_{450\text{nm}}$  (or dual  $\lambda_{450/650\text{nm}}$ , if possible).


## Calculations of Results

Relative optical density (OD) values obtained using CHEMICON's® Transcription Factor Assay Kit may be compared with known standards or other test samples to obtain relative activities. The following graphs represent unstimulated or TNF $\alpha$ -stimulated HeLa nuclear extracts. The data below is for reference use only and the data should not be used to interpret actual assay results.

**Figure 1:** NF $\kappa$ B p50/p65 Transcription Factor Assay. Titration of TNF $\alpha$ -treated HeLa (i.e. 0.2  $\mu$ g/mL for 30 min. at 37°C) nuclear extract and untreated HeLa nuclear extracts. The assay was performed using the SGT510 kit and its various components with normal parameters of the assay as outlined in the assay protocol.



**Figure 2: NFκB p50/p65 Transcription Factor Assay run with 10 μg/well of either untreated or TNFα-treated HeLa nuclear extract (0.2 μg/mL for 30 min. at 37°C) using the assay included controls. The SGT510 assay with its various controls was run under the normal parameters as outlined in the Assay Protocol above. The results show the sensitive specific signal obtained from the assay.**



---

## Troubleshooting

### 1. Weak or no signal in all wells

Possible explanations:

<b><u>Cause</u></b>	<b><u>Solution</u></b>
Omission of a reagent and/or step in the protocol	Make sure to follow the protocol accurately. Perform assay using the positive control extract.
Improper preparation or storage of a reagent	Make sure that all reagents are stored at their proper temperatures. Avoid repeat freeze/thawing of frozen materials.
Reagent Expired	Check the date on reagents to make sure they have not expired.
Insufficient cell extract per well resulting in low sample concentration	Perform cell extract titer to determine optimal concentration. The Enhanced Transcription Factor Buffer is tolerant of variations in sample volume.
Salt concentrations affecting DNA: protein binding	Reduce the amount of extract used in assay or reduce the amount of salt in extraction buffer. Alternatively, perform buffer exchange.
Plate reader not working well	Allow plate reader light to warm up. Check Plate reader.
TMB/E Substrate too cold	Pre-warm the substrate to room temperature or 37°C.

**2. High signal in all wells**

Possible explanations:

<u>Cause</u>	<u>Solution</u>
Improper or inadequate washing of the wells	Follow the protocol as to the times and volume of each wash. Do not use 1x Enhanced Transcription Factor Assay Buffer for the final wash. 1x Transcription Factor Assay Buffer without Blocking Reagent should be used for the final wash.
Improper antibody dilution	Follow the product insert as to the proper antibody dilutions in the 1x Enhanced Transcription Factor Assay Buffer.
Overdeveloping	Decrease the substrate incubation time in the final step of the assay.

**3. High background in sample wells**

**Cause**

Sample concentration too high

Improper antibody dilution

**Solution**

Adjust the concentration by generating a dose response curve with your sample to determine proper concentration of sample to use.

Follow the product insert as to the proper antibody dilutions in the 1x Enhanced Transcription Factor Assay Buffer.

**4. Weak signal in sample wells**

**Cause**

Sample concentration too low

Improper antibody dilution

**Solution**

Not enough cellular extract was used. Increase volume, or prepare fresh sample as weak signal can result from poor sample preparation.

Follow the product insert as to the proper antibody dilutions in the 1x Enhanced Transcription Factor Assay Buffer.

---

## References

1. DiDonato, J.A., Mercurio, F., and Karin, M. (1995). Phosphorylation of I $\kappa$ B $\alpha$  Precedes but is not sufficient for its dissociation from NF- $\kappa$ B. *Mol. Cell. Biol.* **15**:1302-1311.
2. Huang, T and Miyamoto, S (2001). Postrepression activation of NF- $\kappa$ B requires the amino-terminal nuclear export signal of I $\kappa$ B $\alpha$ . *Mol. Cell. Biol.* **21**:4737-4747.
3. Ghosh, S and Karin, M. (2002). Missing Pieces in the NF- $\kappa$ B puzzle. *Cell* **109**:S81-S96.
4. Baeuerle, P.A. and Baltimore, D. (1988). Activation of DNA-binding activity in an apparently cytoplasmic precursor of the NF-kappa B transcription factor. *Cell* **53**:211-2117.
5. Miyamoto, S., Chiao, P.J., and Verma, I.M. (1994). Enhanced I kappa B alpha degradation is responsible for constitutive NF- kappa B activity in mature murine B-cell lines. *Mol. Cell. Biol.* **14**:3276-3282.

---

## **Warranty**

These products are warranted to perform as described in their labeling and in **CHEMICON**<sup>®</sup> literature when used in accordance with their instructions. THERE ARE NO WARRANTIES WHICH EXTEND BEYOND THIS EXPRESSED WARRANTY AND **CHEMICON**<sup>®</sup> DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR WARRANTY OF FITNESS FOR PARTICULAR PURPOSE. **CHEMICON**<sup>®</sup>'s sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of **CHEMICON**<sup>®</sup>, to repair or replace the products. In no event shall **CHEMICON**<sup>®</sup> be liable for any proximate, incidental or consequential damages in connection with the products.

©2003 - 2006: CHEMICON<sup>®</sup> International, Inc. - By CHEMICON<sup>®</sup> International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

**Cat No. SGT510, SGT515, SGT520**

August 2006  
Revision F: 41431