



CE **Authorised Representative: Millipore (UK) Ltd,**
Fleming Road, Kirkton Campus Livingston EH54 7BN UK.
Tel +44 1506 404000
Fax +44 1506 404073

Normal Goat FITC Conjugate

CATALOG NUMBER:	5202i	QUANTITY:	1.0 mL
LOT NUMBER:		CONCENTRATION:	1 mg/mL
IMMUNOGLOBULIN:	Goat IgG FITC Conjugate		
DILUTION:	Optimal dilution or concentration was established using the Polyclonal Rabies Reagent (Cat. No. 5199i) as a positive control.		
FORMAT:	Purified antibody in protein-buffered solution.		
PRESENTATION:	Purified antibody in 0.02M PB, 0.25M NaCl, pH 7.6 contains 0.1% Sodium Azide as a preservative and 15mg/mL BSA as a stabilizer.		
STORAGE/HANDLING:	Maintain at 2°C to 8°C in undiluted aliquots. When stored at 2°C to 8°C, reagent is stable up to the expiration date printed on the label.		

RECOMMENDED

PROTOCOL:

1. Allow the acetone-fixed Positive and Negative Control impression slides to equilibrate to room temperature.
2. If uncoated slides are used, encircle each impression with an ImmuPen or wax pencil.
3. Add the diluted conjugate directly onto control and test slides. Care must be taken to avoid cross-contamination of slides during the staining procedure.
4. Incubate in a humid chamber at 37°C for 30 minutes.
5. Remove excess antibody from the slides by briefly rinsing with PBS. Next, wash the slides for 10 minutes with PBS. Change the PBS once during this time.
6. Rinse briefly with deionized water and allow to air dry.
7. Mount coverslip using mounting fluid containing 10T to 20% glycerol at pH 8.5. (such as Catalog No. 5096)
8. Examine the mounted slides using a fluorescence microscope.

For in vitro diagnostic use. Export only.

Important Note:

During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

©2002 - 2009: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive Temecula, CA 92590 • USA & Canada
+1(800) 437-7500 • Fax: +1 (951) 676-9209
www.millipore.com

Cat. No. 5202i
June 2009
Revision D: 41099i