



PureProteome™ Nickel Magnetic Beads

User Guide

Introduction

The purification of recombinant proteins has been facilitated by the addition of affinity purification tags to the protein sequence. Millipore has developed para-magnetic affinity media for the purification of recombinant, His-tagged proteins based on the well established nickel ion/histidine interaction. Cell lysates prepared under native or denaturing conditions are brought into contact with the PureProteome Nickel Magnetic Beads to bind the His-tagged protein. The beads are isolated using a magnetic rack, followed by wash steps to remove unbound proteins. The bound protein is then eluted at a high purity. This magnetic system is a simple and convenient bench-top format capable of purifying up to a milligram of protein. The PureProteome Nickel Magnetic Beads can be used for screening or for larger scale purification of His-tagged proteins.

Materials Required to Use PureProteome Nickel Magnetic Beads

For optimal performance, the Millipore Magna GriP™ Rack (Cat. no. 20-400) is recommended for use with PureProteome Nickel Magnetic Beads.

Recommended Buffers

Lysis/Equilibration Buffer:

For Purification Under Native Conditions: 50 mM sodium phosphate, 300 mM sodium chloride, 0–10 mM imidazole, pH 8

For Purification Under Denaturing Conditions: 8 M urea, 100 mM sodium phosphate, 10 mM Tris HCl, pH 8

or
6 M Guanidine-HCl, 100 mM sodium phosphate, pH 8

Wash Buffer: 50 mM sodium phosphate, 300 mM sodium chloride, 5–20 mM imidazole (typically 10 mM), pH 8

Elution Buffer: 50 mM sodium phosphate, 300 mM sodium chloride, 300 mM imidazole, pH 8

- NOTES:
1. Lysis/Equilibration and Wash Buffers can be supplemented with 0.1% Tween® 20, Triton® X-100, or Nonidet™ P40 detergents.
 2. When purifying under denaturing conditions, it may be necessary to supplement the Wash and Elution Buffers with either urea or guanidine to prevent protein precipitation.

Sample Preparation Guidelines

This section provides recommendations for sample preparation in order to obtain optimal performance from PureProteome Nickel Magnetic Beads.

Selecting Optimal Conditions

PureProteome Nickel Magnetic Beads may be used for rapid purification of recombinant His-tagged proteins under native or denatured conditions. Optimal conditions for recombinant protein purification depend on the solubility of the protein being studied. Moreover, depending on the secondary structure of the protein and His-tag availability, some proteins bind to the resin better in a native state while others may bind more optimally in a denatured state. While denaturing conditions allow purification of insoluble protein, biological activity will be lost and the protein must be refolded to restore function.

To determine whether a recombinant protein is expressed in a soluble form, perform native lysis. After lysate clarification, analyze the protein content of both the supernatant and the pellet by SDS-PAGE and/or western blotting. If the recombinant protein is insoluble, it should be solubilized in urea or guanidine-HCl buffer prior to purification with the PureProteome Nickel Magnetic Beads. If guanidine is used for protein solubilization, dilute the lysate to 0.1 M guanidine prior to the addition of SDS loading buffer. Alternatively, urea may be substituted for guanidine.

Preparation of Lysates

Cells should be lysed in buffers compatible with the Nickel Magnetic Bead purification protocol. Reagents such as DTT, SDS, and citrate should be avoided or kept to a minimum (see “Compatibility with Commonly Used Purification Reagents” table). Protease inhibitors (1 mM phenylmethanesulfonyl fluoride [PMSF] or commercially available inhibitors compatible with metal affinity purification) are recommended but should not contain > 5 mM EDTA. Cells can be disrupted by addition of lysozyme, sonication on ice, freeze-thaw cycles, using a french press or nebulizer. Since lysozyme tends to bind to metal affinity resins, it may contaminate the final preparation if used in concentrations higher than 0.1 mg/mL.

Clarifying Lysates

Clarify cell lysates prior to protein purification by centrifuging at 10,000 × g for 20–30 minutes to remove genomic DNA, insoluble proteins and cell debris. The addition of Benzonase® nuclease to cell lysates will help reduce viscosity and facilitate downstream purification.

Sample Preparation Guidelines, continued

Optimizing Binding

By using the appropriate amount of beads, His-tagged protein binding can be maximized while minimizing non-specific protein binding. Typical binding capacity is 1–3 µg of recombinant protein per µL of bead suspension. Non-specific binding can be significantly reduced by the addition of low concentrations of imidazole to both the Lysis/Equilibration Buffer (5–10 mM) and Wash Buffer (5–20 mM). Higher concentrations of imidazole may prevent binding of His-tagged proteins. The optimal binding/washing condition for Pure Proteome Nickel Magnetic Beads is at a basic pH (7.5–8), in the presence of 0.3 M NaCl. Refer to the “Recommended Buffers” section for optimized buffer recipes.

Elution

His-tagged proteins can be eluted with 100–500 mM imidazole or low pH. Elution conditions will be dependent on the protein of interest, therefore, optimization will be required.

An eluted protein may be quantified by various protein assay methods. Since most of these methodologies are incompatible with high concentrations of imidazole, removal of the imidazole is necessary. Imidazole removal can be achieved by dilution of the sample, dialysis, or diafiltration. For information on diafiltration, see Millipore’s Ultrafiltration Application and Product Guide.

The addition of low concentrations of nonionic detergents helps to reduce viscosity, solubilize proteins, and accelerate bead collection on the magnet. We recommend supplementing Lysis/Equilibration and Wash Buffers with 0.1% Tween 20, Triton X-100, or Nonidet P40 detergents.

Compatibility with Commonly Used Purification Reagents

The chemicals listed below have no significant impact on protein yield. However, compatibility is protein-dependent. Care should be taken when adding new components to the purification buffers.

| Reagents | | Concentration |
|---------------------|------------------------|-----------------|
| Buffers/salts | MOPS, pH 7.4 | 100 mM |
| | Tris acetate, pH 7.4 | 100 mM |
| | Tris HCl, pH 7.4 | 100 mM |
| | HEPES, pH 7.4 | 100 mM |
| | NaCl | 1 M |
| Reducing reagents | DTT | 5 mM |
| | Beta-mercaptoethanol | 20 mM |
| Chelating reagents | EDTA | up to 5 mM* |
| Detergents | Tween 20 detergent | 2% |
| | Nonidet P40 detergent | 2% |
| | Triton X-100 detergent | 2% |
| | CHAPS | 1% |
| | SDS | Not recommended |
| Denaturing reagents | Urea | 8 M |
| | Guanidine-HCl | 6 M |

* Addition of EDTA may slightly reduce protein binding.

Procedure for Using PureProteome Nickel Magnetic Beads

This protocol is optimized for 10 mL of bacterial culture or 1 mL of prepared lysate. It may be adjusted to account for varying culture/lysate volumes, as well as for protein expression levels.

1. Resuspend magnetic beads by vortexing.
2. Aliquot 200 µL of magnetic bead suspension (200 µL of PureProteome Nickel Magnetic Bead suspension can bind 200–600 µg of His-tagged protein) into a 1.5 mL micro centrifuge tube.
3. Place the tube into the Magna GrIP Rack to collect the beads. Carefully remove storage buffer with a pipette.
4. Resuspend the magnetic beads in 500 µL of Lysis/Equilibration Buffer and incubate with gentle mixing for one minute at room temperature.
5. Place the tube back into the Magna GrIP Rack and remove the buffer.
6. Add 1 mL of cell lysate to the magnetic beads and incubate with gentle mixing for 30 minutes at room temperature.
7. Place the tube back into the Magna GrIP Rack and allow the beads to collect on the magnet. With the tube still seated in the Magna GrIP Rack, invert to remove residual beads from tube cap. Alternatively, residual liquid can be removed from the cap by brief centrifugation. Capture beads and remove lysate.
8. Wash the magnetic beads by incubating in 500 µL Wash Buffer with gentle mixing for one minute at room temperature.
9. Place the tube back into the Magna GrIP Rack and allow the beads to collect on magnet. Remove Wash Buffer.
10. Repeat steps 8 and 9 two more times.
11. Elute the bound protein by adding 100 µL of Elution Buffer. Incubate with gentle mixing for two minutes at room temperature.
12. Place the tube back into the Magna GrIP Rack, allow beads to collect, and transfer the eluted fraction into a clean collection tube.
13. Repeat the elution steps 11 and 12 one more time. The first elution will contain the majority of the recombinant protein. If desired, both elutions can be combined, but this will result in an overall lower protein concentration.

Specifications

| | |
|-----------------|---|
| Matrix | Polymer-coated inorganic beads, pre-charged with Ni ²⁺ |
| Particle form | Spherical |
| Bead diameter | 10 µm (nominal) |
| Storage | 2–8 °C. Do not freeze. |
| Capacity* | Typically 1–3 µg of a 30 kDa His-tagged protein per µL of bead suspension |
| Protein purity* | Typically ≥ 90% by gel densitometry |

* Performance will vary depending on the properties of the individual proteins being used.

Troubleshooting/Optimization

| Problem: Recombinant protein is not present in the eluate | |
|--|---|
| Cause | Solution |
| Protein is insoluble (present in inclusion bodies). | After cell lysate clarification by centrifugation, check both supernatant and pellet for the presence of the recombinant protein. Bind under denaturing conditions. |
| Protein binds to PureProteome Nickel Magnetic Beads tightly and does not elute. | After the binding step, check the lysate for the depletion of the recombinant protein. Increase concentration of imidazole in the elution buffer. |
| Recombinant protein is degraded during the process of cell lysis. | Add protease inhibitors to the cell lysate. Perform sonication in short bursts, keeping the tube on ice between bursts. |
| pH of the Lysis/Equilibration Buffer is incorrect. | Binding buffer pH should be between 7 and 8. Acidic buffer will inhibit binding. |
| SDS or citrate in the Lysis/Equilibration Buffer prevents binding of His-tagged protein. | Refer to "Compatibility with Commonly Used Purification Reagents" section. |
| High concentration of imidazole in the Lysis/Equilibration Buffer. | Do not use more than 5–10 mM imidazole. If the problem persists, use Lysis/Equilibration buffer without imidazole. |
| His-tag is not exposed for binding to the affinity resin. | The protein may require denaturing conditions for binding. Express protein with His-tag on the other terminus. |
| Protein is not expressed correctly or His-tag is not present. | Sequence the plasmid, check cell lysate for the presence of correctly expressed protein via an independent method (e.g., western blotting) |
| Protein induction insufficient. | Confirm that cells have grown to appropriate optical density prior to induction of recombinant protein expression. Make sure appropriate reagents are used for induction of expression. |

Troubleshooting/Optimization, continued

| Problem: Magnetic beads do not collect on the magnet | |
|---|---|
| Cause | Solution |
| Magnet strength is not sufficient. | The Magna GriP Rack is required for optimal performance with PureProteome Nickel Magnetic Beads. |
| Cell lysate is too viscous. | Add nuclease to reduce lysate viscosity. Add 0.1% Tween 20, Triton X-100, or Nonidet P40 detergent to the cell lysate. |
| Problem: High non-specific binding | |
| Cause | Solution |
| Ionic strength of the Lysis/Equilibration and Wash Buffers is not sufficient. | Adjust salt concentration to at least 0.3 M NaCl. Concentrations up to 1 M NaCl are compatible. |
| Binding of non-His tagged proteins. | Add 5-20 mM imidazole to the Lysis/Equilibration and Wash Buffers. Up to 50 mM imidazole may be needed in some cases. |
| Lysozyme contamination of the cell lysate. | Lysozyme tends to bind metal affinity resins. Use lower concentration of lysozyme (0.1 mg/mL), and/or employ other means to disrupt cells (sonication, french press, freeze-thaw cycles, or a nebulizer). |
| Insufficient washing of the beads. | Add wash steps or use a higher volume of Wash Buffer. |
| Magnetic beads carry over. | If some magnetic beads are present in the final elution fractions, non-specifically bound proteins may be stripped off in the SDS loading buffer and appear on the gel. Normally, these proteins will not elute with imidazole. Allow longer time for the beads to collect on the magnet to avoid carry over. |
| His-tagged protein is degraded. | Degradation products also bind PureProteome Nickel Magnetic Beads and may look like contaminating bands on SDS-PAGE gel. Add protease inhibitors to the cell lysis buffer, perform lysis on ice. |

Performance

Purification with PureProteome Nickel Magnetic Beads



Replicate purification of 38 kDa (lanes 2, 3) and 24 kDa (lanes 5, 6) His-tagged proteins using PureProteome Nickel Magnetic Beads. Lanes 1 and 4, corresponding *E. coli* cell lysates; lanes 2, 3, 5, 6, eluted recombinant proteins.

Product Ordering Information

| Description | Qty/Pk | Cat. No. |
|------------------------------------|-----------------------|----------------------------|
| PureProteome Nickel Magnetic Beads | 2 x 1 mL 1 x 10 mL | LSKM AGH 02 LSKM AGH 10 |
| Magna GrIP Rack | 1 | 20-400 |

Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Millipore Corporation or an affiliated corporation. Neither Millipore Corporation nor any of its affiliated corporations assumes responsibility for any errors that may appear in this document. This manual is believed to be complete and accurate at the time of publication. In no event shall Millipore Corporation be liable for incidental or consequential damages in connection with or arising from the use of this manual. It is the purchasers' responsibility to determine whether any other licenses or rights are required in order to perform their specific application.

Technical Assistance

For more information, contact the Millipore office nearest you. In the U.S., call **1-800-MILLIPORE** (1-800-645-5476). Outside the U.S., see your Millipore catalogue for the phone number of the office nearest you or go to our web site at www.millipore.com/ offices for up-to-date worldwide contact information. You can also visit the tech service page on our web site at www.millipore.com/techservice.

Standard Warranty

Millipore Corporation ("Millipore") warrants its products will meet their applicable published specifications when used in accordance with their applicable instructions for a period of one year from shipment of the products. **MILLIPORE MAKES NO OTHER WARRANTY, EXPRESSED OR IMPLIED. THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.** The warranty provided herein and the data, specifications and descriptions of Millipore products appearing in Millipore's published catalogues and product literature may not be altered except by express written agreement signed by an officer of Millipore. Representations, oral or written, which are inconsistent with this warranty or such publications are not authorized and if given, should not be relied upon.

In the event of a breach of the foregoing warranty, Millipore's sole obligation shall be to repair or replace, at its option, the applicable product or part thereof, provided the customer notifies Millipore promptly of any such breach. If after exercising reasonable efforts, Millipore is unable to repair or replace the product or part, then Millipore shall refund to the customer all monies paid for such applicable product or part. **MILLIPORE SHALL NOT BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, SPECIAL OR ANY OTHER INDIRECT DAMAGES RESULTING FROM ECONOMIC LOSS OR PROPERTY DAMAGE SUSTAINED BY ANY CUSTOMER FROM THE USE OF ITS PRODUCTS.**

LSKMAGHXX, Rev. A, 10/08

© 2008. Millipore Corporation. All rights reserved.

Millipore is a registered trademark of Millipore Corporation.
PureProteome, Magna GrIP, and the M mark are trademarks of Millipore Corporation.
Benzonase is a registered trademark of Novagen.
Nonidet is a trademark of Royal Dutch/Shell Group.
Triton is a registered trademark of Union Carbide Corporation.
Tween is a registered trademark of ICI Americas Inc.