



Statement of Animal Origin

Catalog Number: KWSCA04TT3

Item Description: OPTICAP XL 4 W/MILLIGARD-4" 1.2/0.5UM TC/TC 3PK

Based on the information received from our suppliers, some raw materials contain animal derived materials. The information we have received from our suppliers is summarized below:

A raw material contains a tallow-derived material that originates from bovine sourced from the United States and Canada. The supplier has stated the following processing conditions for the tallow-derived material: trans-esterification or hydrolysis at not less than 200 °C for not less than 20 minutes under pressure, and distillation at 200 °C. The processing conditions meet the requirements of the "Note for Guidance on Minimising the Risk of Transmitting Animal Spongiform Encephalopathy Agents via Human and Veterinary Medicinal Products" (EMEA/410/01 Rev. 2).

Fungayi Rowlette

Quality Assurance

Millipore's Quality Compliance includes a Material Sourcing and Animal Origin Program to actively manage its supply chain to ensure compliance with applicable regulatory requirements.

The wording or format of this statement has been modified for clarification only. Please contact Tech Service at 1-800-MILLIPORE with any questions.

Revision: C

Revision date: 03/03/2006

Mfg Location: JF