

**MOUSE ANTI-SV40 T AND t ANTIGEN (N-TERMINUS)
MONOCLONAL ANTIBODY**

CATALOG NUMBER: MAB988

LOT NUMBER:

QUANTITY: 100 µL

SPECIFICITY: Reacts with a determinant on the N-terminal region of the SV-40 large T and small t antigens. Cross reacts with the large T antigen of BK virus (a human papovavirus).

IMMUNOGEN: SV40 transformed, BALB / 3T3 cells.

ISOTYPE: IgG_{2a}

APPLICATIONS: Dilute with buffer pH 7.5 - 8.0 to desired working volumes. For extensive dilution, protein-containing or other stabilizing medium should be used. Final working dilutions must be determined by end user.

FORMAT: Ascites.

STORAGE/HANDLING: Maintain at -20°C for up to 12 months. Avoid repeated freeze/thaw cycles.

REFERENCES: *J. Virology* 57:1168-1172 (1986).

For research use only; not for use as a diagnostic.

© 2002-2004: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (909) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com