

**CHICKEN ANTI-PROTEOLIPID PROTEIN (PLP)
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB15454
LOT NUMBER:	
QUANTITY:	500 μ L
CONCENTRATION:	100 μ g/mL
SPECIFICITY:	PLP [Proteolipid Protein].
IMMUNOGEN:	Synthetic peptide corresponding to human/mouse PLP protein.
APPLICATIONS:	Immunohistochemistry: 1:100-1:200 using 2% paraformaldehyde fixed tissue. Immunocytochemistry: 1:100-1:200 using 2% paraformaldehyde fixed cells. Optimal working dilutions must be determined by end user.
SPECIES REACTIVITIES:	Mouse. Other species have not yet been tested. It is expected that the antibody will also react with human due to immunogen sequence homology.
FORMAT:	Affinity purified chicken immunoglobulin.
PRESENTATION:	Liquid in PBS, pH 7.2 with 0.02% sodium azide.
STORAGE/HANDLING:	Maintain at 2-8°C in undiluted aliquots for up to 6 months after date of receipt.

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000
www.millipore.com