

**RABBIT ANTI-HUMAN CYTOCHROME P450 CYP2B6 and RAT CYP2B1/2
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB1283
LOT NUMBER:	xxxxx
QUANTITY:	100 μ L
SPECIFICITY:	Reacts with cytochrome P450 CYP2B6 in human hepatic microsomal fraction and to CYP2B1 and CYP2B2 in rat hepatic microsomal fraction.
IMMUNOGEN:	Synthetic peptide
APPLICATIONS:	Western blot: CYP2B1 and CYP2B2 can be separated electrophoretically. Immunohistochemistry on frozen and formaldehyde treated sections. Optimal working dilutions must be determined by end user.
FORMAT:	Rabbit antiserum
PRESENTATION:	Liquid, no preservatives.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months after date of receipt. Avoid repeated freeze/thaw cycles.
REFERENCES:	<i>Biochemical Pharmacology</i> (1998) 56 :377-387. <i>Biochemical Pharmacology</i> (1995) 50 :1517-1520.
RELATED REFERENCES:	<i>Federation Proc.</i> (1987) 46 :2563-2566. <i>Adv. Enzymol.</i> (1987) 60 :35-87. <i>DNA and Cell Biology</i> (1993) 12 :1-51.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ l or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

©2003 - 2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com