

**MOUSE ANTI-HUMAN HERPES VIRUS TYPE 6 (HHV6), gp60/110  
MONOCLONAL ANTIBODY**

<b>CATALOG NUMBER:</b>	MAB8537
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	100 µg
<b>CONCENTRATION:</b>	100 µg/mL
<b>SPECIFICITY:</b>	Specific for the gp 60/110 envelope glycoprotein of HHV6.
<b>CELLULAR LOCALIZATION:</b>	Cytoplasmic and membrane
<b>STRAIN RECOGNITION:</b>	Types A & B
<b>ISOTYPE:</b>	IgG <sub>2b</sub>
<b>APPLICATIONS:</b>	Western Blot 1:10 - 1:50 Immunofluorescence 1:10-1:50 Immunoprecipitation 1:10-1:50 Optimal working dilutions must be determined by end user.
<b>FORMAT:</b>	Purified immunoglobulin.
<b>PRESENTATION:</b>	Liquid in 0.01 M phosphate buffered saline (PBS) pH 7.2, containing 0.1% sodium azide.
<b>STORAGE/HANDLING:</b>	Maintain frozen at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

(c) 2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233  
Australia +61 3 9839 2000  
[www.millipore.com](http://www.millipore.com)