

**GABA B RECEPTOR R1  
CONTROL PEPTIDE**

**CATALOG NUMBER:** AG324

**LOT NUMBER:**

**QUANTITY:** 100 µg

**CONCENTRATION:** 1.0 mg/mL

**DESCRIPTION:** Peptide from rat GABA B Receptor R1.

**APPLICATIONS:** For use in blocking the reactivity of AB1531.  
Optimal working dilution must be determined by the end user.

**PRESENTATION:** Liquid in PBS.

**STORAGE/HANDLING:** Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

*For research use only; not for use as a diagnostic.*

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*