



**CHEMICON International, Inc.** 28835 Single Oak Drive • Temecula, CA 92590  
Phone: 909-676-8080 • 800-437-7500 • Fax: 909-676-9209  
custserv@chemicon.com • techserv@chemicon.com • www.chemicon.com

## MOUSE ANTI-SNF2 $\beta$ MONOCLONAL ANTIBODY

**CATALOG NUMBER:** MAB3452

**LOT NUMBER:**

**QUANTITY:** 100  $\mu$ L

**SPECIFICITY:** MAB3452 reacts with SNF2 $\beta$ .

**IMMUNOGEN** Recombinant mouse SNF2 $\beta$  (C-terminal region).

**ISOTYPE:** IgG1

**APPLICATIONS:** Western blot: 1:500-1:5,000  
Immunofluorescence: 1:500-1:5,000  
ELISA: 1:500-1:5,000  
Optimal working dilutions must be determined by the end user.

**SPECIES REACTIVITY:** Mouse. Has not been tested on other species.

**FORMAT:** Ascites

**PRESENTATION:** Liquid. Contains no preservative.

**STORAGE/HANDLING:** Maintain at -20°C in undiluted aliquots up to 12 months. Avoid repeated freeze/thaw cycles.

*For research use only; not for use as a diagnostic.*

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200  $\mu$ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*