

**MOUSE ANTI-HUMAN SURFACTANT PROTEIN A (SP-A)
MONOCLONAL ANTIBODY**

CATALOG NO.: MAB3270

LOT NUMBER:

QUANTITY: 100 μ L

SPECIFICITY: Reacts with human SP-A.

IMMUNOGEN: Human SP-A purified from lavage material.

ISOTYPE: IgM

CLONE: #3

APPLICATIONS: Immunoblot
ELISA
Optimal dilutions must be determined by the end user.

SPECIES REACTIVITY: Human.

FORMAT: Ascites fluid.

PRESENTATION: Liquid. Contains no preservative

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to twelve months. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.