

NORMAL RAT SERUM

CATALOG NUMBER: S24-100 mL

LOT NUMBER: xxxxxxxx

QUANTITY: 100 mL

PRODUCT DESCRIPTION: Normal rat serum collected from healthy normal rats.

FORMAT: Nonsterile liquid. Contains no preservatives.

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

© 2003-2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com