

**RABBIT ANTI- VASP  
POLYCLONAL ANTIBODY**

<b>CATALOG NUMBER:</b>	AB3846
<b>LOT NUMBER:</b>	XXXXXXX
<b>QUANTITY:</b>	100µL
<b>SPECIFICITY:</b>	Recognizes VASP (Vasodilator-stimulated phosphoprotein). VASP belongs to the Ena/VASP family of adaptor proteins linking the cytoskeletal system to the signal transduction pathways and functions in cytoskeletal organization, fibroblast migration, platelet activation and axon guidance.
<b>IMMUNOGEN:</b>	KLH-conjugated, synthetic peptide corresponding to the sequence surrounding Ser157 of human VASP.
<b>HOST:</b>	Rabbit
<b>APPLICATIONS:</b>	Western blotting: 1:1,000 (for best results, incubate membrane with diluted antibody in 5% BSA, 1X TBS, 0.1% Tween-20 at 4°C with gentle shaking, overnight) <i>Optimal working dilutions must be determined by the end user.</i>
<b>SPECIES REACTIVITIES:</b>	Human, mouse, rat, and monkey.
<b>FORMAT:</b>	Protein A and peptide affinity purified immunoglobulin
<b>PRESENTATION:</b>	Liquid in 10 mM sodium HEPES, pH 7.5, 150 mM NaCl, 100 µg/mL BSA, with 50% glycerol.
<b>STORAGE/HANDLING:</b>	Maintain at -20°C for 12 months, do not aliquot.

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

***For research use only; not for use as a diagnostic.***