

**CHICKEN ANTI-VAMP2
AFFINITY PURIFIED
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB5625
LOT NUMBER:	xxxxxx
QUANTITY:	100 µg
CONCENTRATION:	1 mg/mL
SPECIFICITY:	VAMP2 (Vessicle-associated membrane protein 2, synaptobrevin 2).
IMMUNOGEN:	Human VAMP2.
APPLICATIONS:	Western Blot: 1:500-1:2000 Immunocytochemistry: 1:200 on AtT-20 cells. It is expected that the antibody will also work for immunohistochemistry. Specific reactivity has not been confirmed. Optimal working dilutions must be determined by the end user.
SPECIES REACTIVITY:	Human, rat and mouse.
FORMAT:	Affinity purified immunoglobulin.
PRESENTATION:	Liquid in PBS.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
BACKGROUND REFERENCE:	Sudhof, T.C., et al., <i>Nature</i> (1995) 376 .

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

Now Available - Chicken IgY Negative Control Antibody (Cat. No. PP25),
Chicken IgY Blocking Solution (Cat. No. 2180) and Chicken IgY Precipitating Agarose (Cat. No. 2185)

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000
www.millipore.com