

**GALANIN RECEPTOR 3 (GALR3)
CONTROL PEPTIDE**

CATALOG NUMBER: AG314

LOT NUMBER:

QUANTITY: 100 µg

CONCENTRATION: 1.0 mg/mL

DESCRIPTION: Peptide from within the cytoplasmic, C-terminus of mouse GALR3 (1).

APPLICATIONS: For use in blocking the reactivity of AB5650P.
Optimal working dilution must be determined by the end user.

PRESENTATION: Liquid in PBS containing 0.02% sodium merthiolate.

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt.
Avoid repeated freeze/thaw cycles.

RELATED REFERENCE: 1.) Kolakowski LF et al (1998) Gene Accession # AF042783

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*