

**MOUSE ANTI-SNAIL  
MONOCLONAL ANTIBODY**

<b>CATALOG NUMBER:</b>	MAB5494
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	100 $\mu$ L
<b>SPECIFICITY:</b>	Reacts with the snail protein from drosophila.
<b>IMMUNOGEN:</b>	Full length recombinant drosophila snail protein.
<b>ISOTYPE:</b>	IgG <sub>1</sub>
<b>APPLICATIONS:</b>	Western blot: 1:500-1:5,000 Immunocytochemistry: 1:500-1:5,000 ELISA: 1:500-1:5,000 Optimal working dilutions must be determined by end user.
<b>SPECIES REACTIVITY:</b>	Drosophila. Other species have not yet been tested.
<b>FORMAT:</b>	Ascites fluid.
<b>PRESENTATION:</b>	Liquid. Contains no preservative.
<b>STORAGE/HANDLING:</b>	Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
<b>REFERENCES:</b>	Mauhin, V, et al., <i>Nucleic Acids Res.</i> (1993) <b>21</b> :3951-3957.

*For research use only; not for use as a diagnostic.*

*Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200  $\mu$ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*