

**RABBIT ANTI-PROHORMONE CONVERTASE PC1/PC3
POLYCLONAL ANTIBODY**

- CATALOG NUMBER:** AB9722
- LOT NUMBER:**
- QUANTITY:** 100 μ L
- SPECIFICITY:** Reacts with Prohormone Convertase PC1/PC3.
- IMMUNOGEN:** Amino acids 113-158 of human preproconvertase PC1/PC3.
- APPLICATIONS:** Immunohistochemistry/immunocytochemistry: 1:50-1:250
Western blot: 1:500
RIA: 1:8,000
Optimal working dilutions must be determined by end user.
- SPECIES REACTIVITIES:** Human, rat and mouse.
- FORMAT:** Rabbit antiserum, partially purified.
- PRESENTATION:** Liquid, no preservatives.
- STORAGE/HANDLING:** Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
- RELATED REFERENCE:** *J. Biol. Chem.* (1995) **270**:10136-10146.
Mol. Endocrinol. (1996) **10**:331-341.

Important Note: *During shipment, small volumes of antibody will occasionally become entrapped in the seal of the product vial. For antibodies with volumes of 200 μ l or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial on a tabletop centrifuge to dislodge any liquid in the container's cap.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.