
HIV-1 NEF ANTIGEN

CATALOG NUMBER:	AG818	QUANTITY:	100 µg
LOT NUMBER:	xxxxxxx	CONCENTRATION:	1 mg/mL
DESCRIPTION:	HIV-1 NEF antigen (recombinant) a.a 3 to a.a 190 of the HIV-1 nef region 20kDa. Reacts strongly with human I HIV positive serum.		
PURITY:	>95%; Purity of proteins is evaluated by SDS-PAGE, by measuring optical density at 280 nm and by method of Bradford et al.		
PRESENTATION:	8 M urea; 20mM Tris-HCl pH 8.0; 10mM β-mercaptoethanol		
STORAGE/HANDLING:	Store at 2-8°C for up to 6 months.		

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.