

**MOUSE ANTI-VARICELLA ZOSTER
MONOCLONAL ANTIBODY**

CATALOG NUMBER: MAB8616

LOT NUMBER:

QUANTITY: 100 µg

CONCENTRATION: 4.9 mg/mL

SPECIFICITY: Recognizes the VZV immediate early gene 62. No cross-reactivity to HSV.

IMMUNOGEN: VZV HZ strain.

ISOTYPE: IgG₁

APPLICATIONS: Indirect immunofluorescence: shows a nuclear and cytoplasmic staining pattern.
Immunoprecipitation
Western blotting
Final working dilutions must be determined by end user.

FORMAT: Purified immunoglobulin in 0.02M PB with 0.25M NaCl, pH=7.6 with 0.1% Sodium azide.

STORAGE/HANDLING: Maintain at 2-8°C for up to 12 months.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes 200µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

Have you published a paper using a Chemicon primary antibody? You may be eligible for credit toward future purchases of Chemicon products!
Contact our Technical Service department for more information.

For research use only; not for use as a diagnostic.

© 2002-2004: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com