

**SHEEP ANTI-VASOACTIVE INTESTINAL PEPTIDE (VIP)
POLYCLONAL ANTIBODY**

CATALOG NO: AB1581

LOT NUMBER:

QUANTITY: 100 µL

SPECIFICITY: Vasoactive Intestinal Peptide (VIP).

IMMUNOGEN: Synthetic peptide (1-28) coupled to KLH with glutaraldehyde.

APPLICATIONS: Immunohistochemistry: 1:1,000
Optimal working dilutions must be determined by end user.

SPECIES REACTIVITIES: Guinea pig, rat and rabbit.

FORMAT: Sheep serum.

PRESENTATION: Lyophilized, no preservatives. Reconstitute with 100 µL of sterile distilled water. Centrifuge to remove any residue.

STORAGE/HANDLING: Maintain at lyophilized material at -20 to -70°C in undiluted aliquots for up to 12 months after date of receipt. After reconstitution maintain at -20 to -70°C for up to 6 months. Avoid repeated freeze/thaw cycles. Glycerol (ACS grade or better) can be added (1:1) for additional stability.

REFERENCE: Zhu, BS, et al., *J. Comp. Neurol.* (1997) **389**:484-495.

For research use only; not for use as a diagnostic.

© 2003-2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com