

RABBIT ANTI-BETA GALACTOSIDASE POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB986
LOT NUMBER:	
QUANTITY:	2 mL
SPECIFICITY:	Specific for Beta galactosidase from <i>E coli</i> by IEP. May cross react with Beta galactosidase from other species.
IMMUNOGEN:	Beta galactosidase from <i>E-coli</i> .
APPLICATIONS:	Immunoblotting: 1:500-1:2,000 Immunoprecipitation ELISA: 1:65,000 when tested against 1 µg of immunogen by a standard sandwich ELISA. Optimal working dilutions must be determined by the end user.
FORMAT:	Rabbit antisera prepared by delipidation and defibrination.
PRESENTATION:	Lyophilized. Reconstitute with 2 mL of sterile deionized water. Contains 0.02M Potassium Phosphate, 0.15M NaCl, pH 7.2 and 0.01% sodium azide after reconstitution.
STORAGE/HANDLING:	Maintain lyophilized material at 2-8°C for up to 12 months after date of receipt. After reconstitution maintain in undiluted aliquots at -20°C for up to 6 months. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000
www.millipore.com