

CHICKEN ANTI-PYY, N-TERMINUS AFFINITY PURIFIED POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB15668
LOT NUMBER:	
QUANTITY:	100 μ L
CONCENTRATION:	0.5 mg/mL
SPECIFICITY:	PYY, amino acids 3-18.
IMMUNOGEN:	Synthetic peptide from the N-terminus of human PYY (amino acids 3-18).
APPLICATIONS:	ELISA: 1:200-1:1,000 when using as the coating antibody and 1:2,000-1:5,000 when using as the detection antibody in a sandwich ELISA. Optimal working dilutions must be determined by end user.
SPECIES REACTIVITIES:	Human, rat and mouse. Other species have not yet been tested
FORMAT:	Affinity purified chicken IgY.
PRESENTATION:	Liquid in PBS with 0.02% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
IMPORTANT NOTE:	During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.