

HUMAN MYELIN BASIC PROTEIN (MBP)

CATALOG NUMBER:	AG42P
LOT NUMBER:	
QUANTITY:	500 µg
SOURCE:	Spinal cord
PURITY:	Approximately 98% pure as determined by SDS-PAGE
RECONSTITUTION:	Restore lyophilized protein with acidic buffer, pH 4-5. After reconstitution, restore to normal, physiologic pH. Reconstitute only as needed for one week's work. MBPs are very unstable in the reconstituted format.
FORMAT:	Lyophilized, no preservatives.
STORAGE/HANDLING:	Maintain lyophilized material at -20°C for up to 12 months. After reconstitution, store in aliquots at -20°C up to 8 weeks. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.