



CHEMICON International, Inc. 28835 Single Oak Drive • Temecula, CA. 92590
Phone: 909-676-8080 • 800-437-7500 • Fax 909-676-9209
custserv@chemicon.com • techserv@chemicon.com • www.chemicon.com

**MOUSE ANTI-CANINE SURFACTANT PROTEIN A (SP-A)
MONOCLONAL ANTIBODY**

CATALOG NO.: MAB3272

LOT NUMBER:

QUANTITY: 100 μ L

SPECIFICITY: Reacts with canine SP-A.

IMMUNOGEN: Canine SP-A purified from lavage material.

ISOTYPE: Ig

CLONE: 3-4

APPLICATIONS: ELISA
Optimal dilutions must be determined by the end user.

SPECIES REACTIVITY: Canine.

FORMAT: Ascites fluid.

PRESENTATION: Liquid. Contains no preservative

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to twelve months. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.