

**MOUSE ANTI-HUMAN PECAM
MONOCLONAL ANTIBODY**

CATALOG NUMBER: MAB2148

LOT NUMBER:

QUANTITY: 100 µg

CONCENTRATION: 1 mg/mL

SPECIFICITY: Human PECAM

IMMUNOGEN: Human resting Umbilical Vein Endothelial Cells

ISOTYPE: IgG₁

CLONE: P2B1

APPLICATIONS: Immunohistochemistry: for use on acetone fixed tissue
Immunocytochemistry: for use on monocytes (HL60), platelets, endothelial cells, activated T cells and some T cell lines (HUT 78)
Immunoprecipitation
FACS
ELISA
Optimal working dilutions must be determined by end user.

FORMAT: Purified from ascites by protein A chromatography.

PRESENTATION: Buffer of 0.02 M Phosphate, 0.25 M Sodium Chloride, pH 7.6 with 0.1% Sodium Azide.

STORAGE/HANDLING: Store at 2 - 8°C.

Important Note: *During shipment, small volumes of antibody will occasionally become entrapped in the seal of the product vial. For antibodies with volumes of 200 µl or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.