

**GOAT ANTI MOUSE IgG  
AFFINITY PURIFIED  
SECONDARY ANTIBODY**

<b>CATALOG NUMBER:</b>	AP124
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	2 mg
<b>CONCENTRATION:</b>	xx mg/mL
<b>SPECIFICITY:</b>	Reacts with Mouse IgG gamma chain, as well as the light chains from all mouse immunoglobulin classes. Monospecific on IEP against whole mouse serum.
<b>IMMUNOGEN:</b>	Highly purified Mouse IgG isolated from pooled normal mouse serum.
<b>PURITY:</b>	Single band on IEP at a concentration of 10 mg/mL when run against rabbit anti goat whole serum.
<b>FORMAT:</b>	Liquid in 0.02M Phosphate Buffer, 0.25M Sodium Chloride and 0.1% Sodium Azide. Filtered through a 0.45 micron membrane filter.
<b>STORAGE:</b>	Store at +2-8°C in undiluted aliquots for up to six months.

*For research use only; not for use as a diagnostic.*

© 2003-2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233  
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide  
[www.chemicon.com](http://www.chemicon.com) • [custserv@chemicon.com](mailto:custserv@chemicon.com) • [techserv@chemicon.com](mailto:techserv@chemicon.com)