

**MOUSE ANTI-VARICELLA ZOSTER  
MONOCLONAL ANTIBODY**

**CATALOG NUMBER:** MAB8614

**LOT NUMBER:**

**QUANTITY:** 100 µL

**SPECIFICITY:** Recognizes the VZV nucleocapsid by indirect immunofluorescence and immunoprecipitation.  
Shows cross-reactivity to HSV.

**IMMUNOGEN:** VZV HZ strain.

**ISOTYPE:** IgG<sub>1</sub>

**APPLICATIONS:** Indirect immunofluorescence: shows a nuclear staining pattern.  
Immunoprecipitation  
Optimal working dilutions must be determined by the end user.

**FORMAT:** Ascites fluid, with Sodium azide.

**STORAGE/HANDLING:** Maintain at -20°C in aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.

Have you published a paper using a Chemicon primary antibody? You may be eligible for credit toward future purchases of Chemicon products!  
Contact our Technical Service department for more information.

*For research use only; not for use as a diagnostic.*

© 2002-2004: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233  
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide  
[www.chemicon.com](http://www.chemicon.com) • [custserv@chemicon.com](mailto:custserv@chemicon.com) • [techserv@chemicon.com](mailto:techserv@chemicon.com)