

**HUMAN PROSTATE SPECIFIC ANTIGEN (PSA)
PURIFIED PROTEIN**

CATALOG NUMBER: AG650

LOT NUMBER:

QUANTITY: 100 µg

CONCENTRATION: 1 mg/mL

SOURCE: Human seminal fluid. Sera from all donors was tested and found to be negative for syphilis, HBsAg, anti-HBs, anti-HCV and antibodies against HIV 1 and 2. No FDA approved test exists for seminal fluid at this time.

PURITY: 98% or greater, determined by SDS-PAGE.

FORMAT/STORAGE: Liquid in PBS, no preservative.

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

REFERENCE: Zisman, A. et. al., *J. Urology* (1995) **154**(3): 1052-1055.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.