



**ANTI-MOUSE IgG, H&L chains,
AFFINITY PURIFIED, FITC CONJUGATED,
POLYCLONAL ANTIBODY
CERTIFICATE OF ANALYSIS**

CATALOG NUMBER:	AP326F	QUANTITY:	2 mL
LOT NUMBER:		CONCENTRATION:	
HOST/ISOTYPE:	Sheep, IgG		
SPECIFICITY:	Prior to conjugation, the antibody was specific for mouse IgG by immunoelectrophoresis against mouse serum. Undersired immune reactivity has been removed by solid phase adsorption against mouse IgM and IgA. Cross-reactivity to human and rabbit serum proteins is determined by ELISA and found to be less than 0.05%		
IMMUNOGEN:	Purified mouse IgG immunoglobulins from pooled mouse serum.		
APPLICATIONS:	IH:1:10-1:400 ELISA: 1:100-1:1000 <i>Optimal dilutions to be determined by enduser and exact assay.</i>		
SPECIES REACTIVITY:	Mouse, other species not examined.		
PRESENTATION:	Liquid in 20 mM Tris, 70 mM NaCl pH 8.0 with 0.2% BSA and 0.1% sodium azide as a preservative.		
STORAGE/HANDLING:	Store at 2-8°C protected from light and temperature fluxuations for up to the expiration date given on the product label.		

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals
©2002 - 2009: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

28820 Single Oak Drive • Temecula, CA 92590
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com