

**RABBIT ANTI-VAV1
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB4094
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	0.58 mg/mL
SPECIFICITY:	Recognizes Vav-1. The Vav-1 proto-oncogene is a member of the Dbl family of guanine nucleotide exchange factors (GEFs) for the Rho family of small GTPases. The protein is important in hematopoiesis, playing a role in T-cell and B-cell development and activation.
IMMUNOGEN:	Synthetic peptide corresponding to amino acids near the C-terminus of human Vav-1.
APPLICATIONS:	Western blot: 1.0 µg/mL Optimal working dilutions must be determined by end user.
SPECIES REACTIVITY:	Human. Reactivity with other species has not been determined.
FORMAT:	Affinity purified immunoglobulin
PRESENTATION:	Liquid in 0.02M phosphate buffer, 0.25M NaCl, with 0.1% sodium azide as a preservative.
STORAGE/HANDLING:	Maintain at 2° to 8°C in undiluted aliquots for up to one year from date of receipt.

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

© 2004: CHEMICON[®] International, Inc. - By CHEMICON[®] International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.