

**RABBIT ANTI-SERINE/THREONINE PROTEIN
PHOSPHATASE 6C (V/C) POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB4091	QUANTITY:	100 µg
LOT NUMBER:		CONCENTRATION:	1 mg/ml
		HOST/ISOTYPE:	Rabbit IgG
SPECIFICITY:	Recognizes the catalytic subunit of Protein Phosphatase 6, MW ~33kDa.		
IMMUNOGEN:	The immunogen used to the purified peptide conjugated to KLH corresponding to the sequence NH ₂ -Cys-Ala-Val-Pro-Asp-Ser-Glu-Arg-Val-Ile- Pro-Pro-Arg-Thr-Thr-Thr-Pro-Tyr-COOH.		
APPLICATIONS:	Immunoblotting (1-10 ug/ml) Immunoprecipitation (10-20ug/ml) Optimal working dilutions must be determined by end user.		
SPECIES REACTIVITY:	Human, Mouse, Rat, Bovine.		
FORMAT:	Purified immunoglobulin.		
PRESENTATION:	Provided as 0.2µm sterile filtered solution in phosphate buffered saline with 0.08% sodium azide.		
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months from date of receipt. Avoid repeated freeze/thaw cycles.		

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002-2008: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000
www.millipore.com

7-26-04/AB4091/MH/03/24/08/ft