

## RABBIT ANTI-NEURTURIN POLYCLONAL ANTIBODY

<b>CATALOG NUMBER:</b>	AB5404
<b>LOT NUMBER:</b>	xxxxx
<b>QUANTITY:</b>	50 µL
<b>SPECIFICITY:</b>	Neurturin. No cross reactivity with GDNF.
<b>IMMUNOGEN:</b>	Recombinant human Neurturin.
<b>APPLICATIONS:</b>	Immunohistochemistry: 1:1,000-1:2,000 Optimal working dilutions must be determined by end user.
<b>SPECIES REACTIVITIES:</b>	Human, mouse, rat. Other species have not yet been tested.
<b>FORMAT:</b>	Rabbit serum
<b>PRESENTATION:</b>	Lyophilized. Contains no preservative. Reconstitute with 50 µL of sterile distilled water. Centrifuge if necessary to remove any residue.
<b>STORAGE/HANDLING:</b>	Maintain lyophilized material at -20°C to -70°C for up to 6 months. After reconstitution maintain frozen at -20°C in undiluted aliquots for up to 6 months. Avoid repeated freeze/thaw cycles.
<b>RELATED REFERENCE:</b>	Kotzbauer, P.T. et al., (1996) <i>Nature</i> . <b>384</b> :467-470.

*For research use only; not for use as a diagnostic.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2008: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233  
Australia +61 3 9839 2000  
[www.millipore.com](http://www.millipore.com)