



**CHEMICON International, Inc.** 28835 Single Oak Drive • Temecula, CA 92590  
Phone: 909-676-8080 • 800-437-7500 • Fax: 909-676-9209  
custserv@chemicon.com • techserv@chemicon.com • www.chemicon.com

## ANTI-HCV NS4 MONOCLONAL ANTIBODY

**CATALOG NUMBER:** MAB8683

**LOT NUMBER:** 18040241

**QUANTITY:** 100 µg

**SPECIFICITY:** HCV NS4 antigen

**IMMUNOGEN:** Produced against recombinant protein

**CLONE NAME:** 1G7

**APPLICATIONS:** ELISA: 1:1000  
Optimal working dilutions must be determined by end user.

**FORMAT:** Affinity purified

**PRESENTATION:** Liquid at 1 mg/mL in phosphate buffered saline, pH 7.2 with 0.01% sodium azide.

**STORAGE:** Store at -20°C in undiluted aliquots. May be stored at 2-8°C for short term use. Avoid repeated freeze-thaw cycles.

*For research use only; not for use as a diagnostic.*

Chemicon offers a comprehensive line of secondary antibodies, including fluorescein, rhodamine, peroxidase, alkaline phosphatase, Cy3 and Cy5 conjugates. We also offer species-absorbed conjugates for dual labeling to minimize nonspecific background staining. Ask your Chemicon Customer Service Representative for more information about Chemicon's secondary antibody conjugates today!