

ALPHA-SYNUCLEIN, DELTA NAC (Recombinant, *E. coli*)

CATALOG NUMBER: AG551

LOT NUMBER:

QUANTITY: 500 µg

DESCRIPTION: Alpha-synuclein, delta NAC. Alpha-synuclein that lacks the NAC region (amino acids 61-95). Six amino acids were added as a linker.

SEQUENCE: MDVFMKGLSKAKEGVVAAAEEKTKQGVAEAAGKTKEGVLYVGSKTKEGVVHG
VATVAEKTGTEIWMKKDQLGKNEEGAPQEGILEDMPVDPDNEAYEMPSEEGYQDYEP

MOL. WT.: 11,935

PURITY: >95% by SDS-PAGE.

COUNTER ION: Lyophilized in final buffer of 20 mM Tris/HCl, pH 7.5, 100 mM NaCl.

PROTEIN CONTENT: 500 µg per vial

PRESENTATION: White lyophilized powder. Resuspend in sterile water at concentration of 1 mg/ml. This will give you a final of 20 mM Tris/HCl, pH 7.5, 100 mM NaCl.

STORAGE/HANDLING: Maintain lyophilized material at -20°C for up to 12 months after date of receipt. After reconstitution maintain at -20°C to -70°C for up to 2 weeks in undiluted aliquots. Avoid freeze/thaw cycles to avoid aggregation.

For research use only; not for use as a diagnostic.

© 2007: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com