

**MOUSE ANTI-HUMAN SURFACTANT PROTEIN B (SP-B)
MONOCLONAL ANTIBODY**

CATALOG NO.: MAB3276

LOT NUMBER: xxxxxxxx

QUANTITY: 100 μ L

SPECIFICITY: Reacts with human SP-B.

IMMUNOGEN: Human SP-B.

ISOTYPE: IgG_{2b}

CLONE: 1B9

APPLICATIONS: Immunoblot under non reducing conditions.
ELISA
Optimal dilutions must be determined by the end user.

SPECIES REACTIVITY: Human.

FORMAT: Ascites fluid.

PRESENTATION: Liquid. Contains no preservative

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to twelve months. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

© 2002-2004: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (909) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com