



CHEMICON International, Inc. 28835 Single Oak Drive • Temecula, CA 92590
Phone: 909-676-8080 • 800-437-7500 • Fax: 909-676-9209
custserv@chemicon.com • techserv@chemicon.com • www.chemicon.com

HIV - 1 GAG ANTIGEN

CATALOG NUMBER: AG817

LOT NUMBER: 19040304

QUANTITY: 100 µg

CONCENTRATION: 1 mg/mL

DESCRIPTION: HIV-1 gag antigen (recombinant) (a.a. 77-436) of the HIV-1 gag region 39kDa. Reacts strongly with human I HIV positive serum.

PURITY: >95%; Purity of proteins is evaluated by SDS-PAGE by measuring optical density at 280 nm and by method of Bradford et al.
Origin: *E. coli*, recombinant
Fusions: β-galactosidase (114kDa) fused at the N-terminus

FORMAT: 8M urea, 20mM Tris-HCl, pH 8.0, 10mM β-mercaptoethanol.

STORAGE/HANDLING: Store at 2-8°C for up to 6 months.

For research use only; not for use as a diagnostic.

Chemicon offers new Tube-O-DIALYZER sets for dialysis of samples as small as one microliter, with 100% sample recovery.
For more information about Tube-O-DIALYZER sets, contact your Chemicon Customer Service Representative today!