

**RABBIT ANTI-CHLORIDE CHANNEL, CLC-3  
AFFINITY PURIFIED  
POLYCLONAL ANTIBODY**

<b>CATALOG NUMBER:</b>	AB5162-200UL
<b>LOT NUMBER:</b>	
<b>QUANTITY:</b>	200 µL
<b>CONCENTRATION:</b>	0.8 mg/mL (after reconstitution)
<b>SPECIFICITY:</b>	Recognizes a full-length CLC-3 voltage-gated Cl protein. The epitope exhibits partial homology with CLC-4 and CLC-5 proteins but not with other known proteins tested so far.
<b>IMMUNOGEN:</b>	GST fusion protein and a portion of rat CLC-3 protein (amino acids 592-661) (SwissProt Accession P51792).
<b>APPLICATIONS:</b>	Western blot: 1:200 using ECL on rat brain membranes. Immunohistochemistry Immunocytochemistry Inhibition of channel activity Dilutions should be made using a carrier protein such as BSA (1-3%) Optimal working dilutions must be determined by the end user.
<b>CONTROL ANTIGEN:</b>	Included free of charge with the antibody is 120 µg of control antigen (lyophilized powder). The stock solution of the antigen can be made up using 100 µL of PBS. For positive control, in Western blot using 20 ng of protein per minigel lane. For negative control, preincubate 3 µg of fusion protein with 1 µg of antibody for one hour at room temperature. Optimal concentrations must be determined by the end user.
<b>SPECIES REACTIVITIES:</b>	Rat. It is expected that the antibody will also react with mouse (70/70), human, rabbit and guinea pig (69/70), <i>Xenopus laevis</i> (61/70) due to sequence homology, Other species have not been tested.
<b>FORMAT:</b>	Affinity purified immunoglobulin.
<b>PRESENTATION:</b>	Lyophilized from phosphate buffered saline, pH 7.4, containing 1% BSA, 5% sucrose as a stabilizer and 0.025% sodium azide as a preservative. Reconstitute with 200 µL of sterile deionized water. Centrifuge antibody preparation before use (10,000 xg for 5 min).
<b>STORAGE/HANDLING:</b>	Maintain lyophilized material at -20°C for up to 12 months. After reconstitution maintain at -20°C in undiluted aliquots for up to 6 months. Avoid repeated freeze/thaw cycles.

**SUGGESTED WESTERN BLOT PROTOCOL**

1. Mix the samples with sample-buffer X 2, and heat 10 min at 70°C.
2. 5-50 µL applied to Minigel lane (0.75-1.5 mm width) and run at standard conditions. (60 mA for 2 1.5 mm Minigel gels, 1.4 h)
3. Transfer in semi-dry system under standard conditions (3 h 200 mA for two minigel gels)
4. Stain the transferred bands with Chemicon BLOT-*FastStain* (Catalog Number 2076).
5. Destain with deionized water.
6. Block with 5% non-fat milk (Marvel or Carnation) in PBS, and 0.025 % sodium azide, overnight at 4°C. The non-fat milk should be dissolved freshly, centrifuged 10,000 rpm for 10 min, and filtered through glass filter (Gelman Acrodisc).
7. Incubation with first antibody 2 h at room temperature or overnight at 4°C in blocking solution. The antibody preparation should be centrifuged before use (10,000 g 5 min.). Optimal working dilutions and incubation time will need to be determined by the end user.
8. Wash 4 x 15 min. with PBS-0.1% tween 20. From this stage, azide should be omitted.
9. Incubation with the secondary antibody (HRP-conjugated goat anti-rabbit antibody, for example Chemicon Catalog Number AP132P, diluted appropriately) 1 h at room temperature.
10. Wash 4 x 15 min. with PBS-0.1% tween 20.
11. Perform ECL with commercial kits (*ChemLUCENT*, Chemicon Catalog Number 2600).



Western blotting of rat brain membranes.

1. AB5162 antibody (1:200)
2. AB5162, preincubated with a control antigen.

*For research use only; not for use as a diagnostic.*

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233  
Australia +61 3 9839 2000  
www.millipore.com

5-10-2002/AB5162-200UL/CB  
/11/02/07/ft