

**RABBIT ANTI-V5
AFFINITY PURIFIED
POLYCLONAL ANTIBODY**

CATALOG NUMBER: AB3792

LOT NUMBER:

QUANTITY: 100 µg

CONCENTRATION: 1 mg/mL

SPECIFICITY: V5 (GKPIPPELLGLDST).

IMMUNOGEN: GKPIPPELLGLDST (V5)

APPLICATIONS: Western blot: 1:1,000-1:10,000 for colorimetric detection; 1:1,000-1:30,000 for chemiluminescent detection.
Immunocytochemistry: 1:100-1:400.
ELISA: 1:100-1:500 for coating; 1:1,000-1:30,000 for primary.
Optimal working dilutions must be determined by the end user.

FORMAT: Affinity purified immunoglobulin

PRESENTATION: Liquid in PBS, pH 7.2 with 0.1% sodium azide.

STORAGE/HANDLING: Maintain at 2-8°C in undiluted for up to 6 months from date of receipt.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

© 2002-2004: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com