

**MOUSE ANTI-SNF2h  
MONOCLONAL ANTIBODY**

**CATALOG NUMBER:** MAB3644

**LOT NUMBER:**

**QUANTITY:** 100  $\mu$ L

**SPECIFICITY:** Recognizes SNF2h (Smarca5).

**IMMUNOGEN:** Synthetic peptide from amino acids 360-374 of human SNF2h.

**ISOTYPE:** IgG<sub>1</sub>

**APPLICATIONS:** Western blot: 1:500-1:5,000  
Immunocytochemistry: 1:500-1:5,000  
ELISA: 1:500-1:5,000  
Optimal working dilutions must be determined by the end user.

**SPECIES REACTIVITIES:** Human and mouse. Reactivity with other species has not been confirmed.

**FORMAT:** Ascites.

**PRESENTATION:** Liquid. Contains no preservative.

**STORAGE/HANDLING:** Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

*For research use only; not for use as a diagnostic.*

**Important Note:** *During shipment, small volumes of antibody will occasionally become entrapped in the seal of the product vial. For antibodies with volumes of 200  $\mu$ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*