

**RABBIT ANTI-DYNAMITIN (p50), aa 1-83
AFFINITY PURIFIED
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB5869P
LOT NUMBER:	xxxxxx
QUANTITY:	50 μ L
CONCENTRATION:	15 μ g/150 μ L
SPECIFICITY:	Dynamitin (p50), aa 1-83
IMMUNOGEN:	Synthetic peptide from human Dynamitin, aa 1-83
APPLICATIONS:	Western blot: 1:1,000 on mouse brain lysate. It is expected that the antibody will also work for immunocytochemistry and immunohistochemistry. Suggested starting 1:50. Optimal working dilutions must be determined by the end user.
SPECIES REACTIVITY:	Human, rat, mouse and rabbit. Other species have not been tested.
FORMAT:	Affinity purified immunoglobulin.
PRESENTATION:	Liquid in PBS with 50% Glycerol, pH 7.3 and 0.1% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in undiluted for up to 6 months from date of receipt. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002-2005: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.