

**RABBIT ANTI-SYNAPTOPHYSIN
AFFINITY PURIFIED
POLYCLONAL ANTIBODY**

CATALOG NUMBER:	AB9272		
LOT NUMBER:			
QUANTITY:	500 μ L		
SPECIFICITY:	Synaptophysin		
IMMUNOGEN:	Synthetic peptide from human synaptophysin.		
APPLICATIONS:	Western blot Immunohistochemistry on formalin fixed, paraffin embedded tissue: 1:300 for 10 minutes at room temperature. Staining of formalin-fixed tissues requires boiling of tissues in 10 mM citrate buffer, pH 6.0 for 10 minutes followed by cooling at room temperature for 20 minutes. Optimal working dilutions must be determined by the end user.		
POSITIVE CONTROL:	Pancreas or Pheochromocytoma	STAINING PATTERN:	Cytoplasmic
SPECIES REACTIVITY:	Human and rat. Other species have not been tested.		
FORMAT:	Affinity purified immunoglobulin.		
PRESENTATION:	Liquid in 10 mM PBS, pH 7.4 containing BSA and 0.1% sodium azide.		
STORAGE/HANDLING:	Maintain at 2-8°C for up to 6 months after date of receipt.		

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 μ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.