

RABBIT ANTI-SUV39H1 POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB3353
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	0.25 mg/mL
SPECIFICITY:	Reacts with SUV39H1.
IMMUNOGEN:	GST fusion protein from the N-terminal of human SUV39H1.
APPLICATIONS:	Western Blot: 1:100-1:500. Predicted molecular weight 47 kDa. ELISA (direct): 1:1,000 Optimal working dilutions must be determined by the end user.
SPECIES REACTIVITY:	Human. Other species have not been tested.
FORMAT:	Purified immunoglobulin.
PRESENTATION:	Liquid in PBS containing 0.09% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
RELATED REFERENCES:	<ol style="list-style-type: none">1. Feng Q, et al. <i>Curr Biol.</i> (2002) 12:1052.2. Huang S. <i>Nat. Rev. Cancer.</i> (2002) 2:469-4763. Berger SL. <i>Curr Opin Genet Dev.</i> (2002) 12:142.4. Wang H, et al. <i>Mol Cell.</i> (2001) 8:1207.5. Ben-Porath I and Cedar H. <i>Mol Cell.</i> (2001) 8:933.6. Du Y, et al. <i>Cancer Res.</i> (2001) 61:8094.7. Peters AH, et al. <i>Cell.</i> (2001) 107:323.8. Vandel L, et al. <i>Mol Cell Biol.</i> (2001) 21:6484.

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*