

RABBIT ANTI-HUMAN VASCULAR ENDOTHELIAL GROWTH FACTOR (VEGF) POLYCLONAL ANTIBODY

CATALOG NUMBER:	CBL42	LOT NUMBER:	
QUANTITY:	500 µg		
ALTERNATE NAMES:	VEGF		
APPLICATIONS:	This antibody is suitable for studies in western blotting, ELISA based assays, RIA and immunoprecipitation. Optimal working dilutions must be determined by the end user.		
SPECIES REACTIVITY:	Reacts with Human. Reactivity with other species has not been determined.		
IMMUNOGEN:	MAP peptide of the N- terminal sequence (aa 1-20) of human VEGF		
FORMAT:	Purified by affinity chromatography using the immunising sequence		
PRESENTATION:	The antiserum is presented as 0.5mg containing 10mM sodium azide. We recommend that each laboratory determine an optimum working titre for use in its particular application.		
STORAGE/HANDLING:	For use within 1 month of purchase store at +4°C, for long term storage aliquot antibody into small volumes and store at -20°C.		

For research use only; not for use as a diagnostic

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

©2002-2003:CHEMICON® International, Inc. - By CHEMICON® International, Inc. All rights reserved. No part of these works may be reproduced in any form without permission in writing.