

**NMDAR1
CONTROL PEPTIDE**

CATALOG NUMBER: AG344

LOT NUMBER:

QUANTITY: 100 µg

CONCENTRATION: 1 mg/mL

DESCRIPTION: Synthetic peptide corresponding to the C-terminus of rat NMDA receptor 1 subunit.

APPLICATIONS: For use in blocking the reactivity of AB1516.
Optimal working dilution must be determined by the end user.

PRESENTATION: Liquid in phosphate buffered saline. Containing no preservative.

STORAGE/HANDLING: Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

©2003 – CHEMICON[®] International, Inc. - By CHEMICON[®] International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.