



**CHEMICON International, Inc.** 28835 Single Oak Drive • Temecula, CA. 92590  
Phone: 909-676-8080 • 800-437-7500 • Fax 909-676-9209  
custserv@chemicon.com • techserv@chemicon.com • www.chemicon.com

**GOAT ANTI-HUMAN SURFACTANT PROTEIN A (SP-A)  
POLYCLONAL ANTIBODY**

**CATALOG NO.:** AB3422

**LOT NUMBER:**

**QUANTITY:** 100  $\mu$ L

**SPECIFICITY:** Reacts strongly with human SP-A.

**IMMUNOGEN:** Human SP-A from human alveolar proteinosis lavage fluid.

**APPLICATIONS:** Immunoblot: 1:1,000-1:5,000  
ELISA: 1:500-1:1,000  
Optimal dilutions must be determined by the end user.

**SPECIES REACTIVITY:** Human

**FORMAT:** Purified immunoglobulin

**PRESENTATION:** Liquid in PBS. Contains no preservative.

**STORAGE/HANDLING:** Maintain at  $-20^{\circ}\text{C}$  in undiluted aliquots for up to twelve months. Avoid repeated freeze/thaw cycles.

**REFERENCE:** *J. Clin. Invest.* (1991) **88**:1976-1981.

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200  $\mu$ L or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

*For research use only; not for use as a diagnostic.*