

**GOAT ANTI-VACCINIA  
POLYCLONAL ANTIBODY**

<b>CATALOG NUMBER:</b>	AB8303
<b>LOT NUMBER:</b>	xxxxxxx
<b>QUANTITY:</b>	100 µg
<b>CONCENTRATION:</b>	23.9 mg/mL
<b>SPECIFICITY:</b>	Reactive to vaccinia virus, variola major, and other pox strains of virus.
<b>IMMUNOGEN:</b>	Vaccinia virus, Lister strain.
<b>APPLICATIONS:</b>	ELISA  Optimal working dilutions must be determined by end user.
<b>FORMAT:</b>	Purified immunoglobulin.
<b>PRESENTATION:</b>	Liquid in PBS with 0.1% sodium azide. Spin briefly prior to use to ensure complete recovery of purified immunoglobulin.
<b>STORAGE/HANDLING:</b>	Maintain at 2-8°C in undiluted aliquots for up to 6 months after date of receipt.

**Important Note:** *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

*For research use only; not for use as a diagnostic.*

© 2003-2005: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (951) 676-9209 • Europe +44 (0) 23 8026 2233  
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide  
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com