

**MOUSE ANTI-DENGUE COMPLEX
MONOCLONAL ANTIBODY**

PART NUMBER:	MAB8705
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	1 mg/mL
SPECIFICITY:	Reacts with all members of the Dengue complex.
IMMUNOGEN:	Dengue type 4 virus antigens.
ISOTYPE:	IgG _{2a}
APPLICATIONS:	Tested for IF and IC. Suggested working dilution: 1:200-400. Final working dilutions must be determined by end user.
FORMAT:	Purified Immunoglobulin in 0.02M PB, 0.25M NaCl, 0.1% sodium azide pH 7.6.
STORAGE/HANDLING:	Maintain at 2-8°C
REFERENCES:	<i>PNAS</i> (2007) 104: 9422-9427. <i>Am. J. Trop. Med. Hyg.</i> (1982) 31:830-836.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µl or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.