

GALANIN CONTROL PEPTIDE

CATALOG NUMBER:	AG370
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	1 mg/mL
DESCRIPTION:	Galanin control peptide.
APPLICATIONS:	Antibody blocking (AB5909) ELISA Optimal working concentration must be determined by the end user.
PRESENTATION:	Liquid in PBS, pH 7.5 containing 0.05% sodium azide.
STORAGE/HANDLING:	Store at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µl or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

© 2002-2004: CHEMICON International, Inc. - By CHEMICON International, Inc. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.

USA & Canada • Phone: +1(800) 437-7500 • Fax: +1 (909) 676-9209 • Europe +44 (0) 23 8026 2233
Australia +61 3 9839 2000 • Germany +49-6192-207300 • ISO Registered worldwide
www.chemicon.com • custserv@chemicon.com • techserv@chemicon.com